

KIT 8080 NYLON TPU ALTERNATING PRESSURE MATTRESS + PUMP



Product Code	28582
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398278
CND	V080701
EAN/UPC	8023279285826

Description

KIT 8080 ALTERNATING PRESSURE MATTRESS + PUMP

INTENSIVE CARE STAGE I, II, III, IV

Advanced systems of alternating pressure mattresses.

The multi-mode pump offers:

- alternating mode: the pressure is periodically redistributed by inflating and deflating cells in order to mimic natural sleep movement patterns and avoid long term pressurization of tissue.
- constant static mode: the pressure is about 2/3 lowered when compared to the same pressure level from alternating mode by means of immersion and envelopment. Visual alarms notify the caregiver of unusual situation, reducing the risk of human error. Suitable for the prevention and treatment of all stages of pressure ulcers.

Manual in: GB, FR, IT, ES, PT, DE, PL, RO, SE, GR.

Cotton lining with zipper

Professional use

CPR straps quickly deflate the mattress for emergencies

ALTERNATING FUNCTION

The pump inflates and deflates the cells with cyclicity, to change the contact position between the patient and the bed, to promote the blood circulation, to prevent and recover the bedsore, reducing the patient's pain and the nursing labor.

Pillow function: the cells of the head position don't alternate, for more comfort to the patient.

Technical Specifications

PUMP:

Size: 30x13.5x11 cm

Cycle time: 12 min

Operating voltage: 220-240 V, 50 Hz

Working pressure: 55-120 mmHg

Case material: Flame retardant ABS

Air output: 7-8 liter/min

TUV-CE IEC/EN60601-1; 60601-1-2

Quick connector with transfer cap to transform mattress into a sealed static surface to facilitate patient transportation.

MATTRESS:

Total size: 200x86x21 cm

Alternating 1 in 2 (group A and B)

Cell height: 20 cm (8")
Cover material: Nylon+PU with cotton lining
Thickness: 0.25 mm
Number of cells: 20+1 spare
Cells material: nylon TPU
Load: 200 kg