

CRYOCLEAR CRYOSURGICAL DEVICE with 16 g cartridge



Product Code	30587
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2450883
CND	M9001
EAN/UPC	0850077006487
GMDN	44724

Description

CRYOCLEAR® WITH 16 G CARTRIDGE - non-refillable

CryoClear® uses Carbon Dioxide which freezes at -79°C to a depth of no more than 2 mm into the skin, delivering a safe and clinically effective treatment to superficial skin indications.

Over 95% high efficacy combined with less destruction of healthy tissue makes it perfect for treating age spots, sun spots, and skin tags.

It offers a better and faster cryotherapy solution, even as combination treatment to chemical peels, facials, micro-needling or dermaplaining.

No contact device for use on different patients at 0.5 cm distance.

Up to 300 seconds treatment time per device (approx 60 treatments of 5 seconds).

Suggested for use at aesthetics and cosmetic sites.

Multilingual pack and instructions (GB, FR, IT, ES, DE, NL, PL, FI).

Made in USA.

Technical Specifications

Freeze time range (sec) - Suggested freezing time as reported in literature:

Actinic Keratosis (sun spots) 5-10

Lentigo (Age Spots) 5-10

Skin Tags 5-10