



Product Code	33680
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2824208
CND	Z12150102
EAN/UPC	8023279336801
PARAF	986897611

## Description

Peak flow meter - 50-400 l/min - pediatric

Regular use of a Peak flow meter can help keep tabs on your asthma by detecting airway narrowing even before you feel any symptoms, giving you time to adjust your medication before your symptoms get worse. Suitable for children as young as 5.

Main features:

- measure day-to-day changes in your breathing
- track the control of your asthma over time
- show how well your treatment is working
- recognize signs of a flare-up before symptoms appear
- know what steps to take when you have signs of an asthma flare-up
- decide when to call your doctor or get emergency care

Gima Peak flow meter can be used together with "Mobile PEF" App.  
Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, SE.

Green area: (normal) PEF >80%

Yellow area: (warning) PEF =60%, <80%

Red area: (danger) PEF <60%