

GERMOCID PERACETIC DISINFECTANT POWDER - bag 16 g



Product Code	36639
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	4266
CND	D99
EAN/UPC	8023279366396
PARAF	977815846

Description

GERMOCID PERACETIC POWDER - bag 16 g - box of 50 bags

Cold sterilizer disinfectant: the peracetic acid which forms from the reaction of sodium percarbonate and tetracetylene diammine (obtained after dissolving the product into water) carries on an efficacious biocidal action against spores, virus (included HIV, HBV, HCV), bacteria, mycobacteria and fungi. At the concentration of 1.6% of product, it carries on sterilization in 10 minutes.

Active in 15 minutes from dissolution and remains active for max 24 hours. Suitable for a cold rapid sterilization and disinfection of endoscopes, surgical and dental instruments, transducer, etc.

Also active in presence of organic materials (previously clean the device with suitable detergent, better with enzymatic base).

Suitable for automatic washing machines. No environmental impact: once in contact with discharge water, it degrades immediately into acetic acid, water, oxygen.

Label/instruction: GB, FR, IT.

Made in Italy.

Technical Specifications

Chemical composition:

Peracetic acid precursors and coformulants.