



Product Code	36815
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	2760930
CND	D99
GTIN	18009110026664

Description

BARRYCIDAL "30 PLUS" DILUTED 2% - bottle 1 l - box of 12

Ready-to-use solution for high-level disinfection of invasive medical devices and surgical instruments. It is used for the treatment of thermolabile and non-autoclavable medical equipment (e.g. ultrasonic transducers, pulse oximeter probes, stethoscopic and endoscopic probes, ECG devices, incubators, equipment for blood pressure measure, radiologic instruments, monitors etc). Quick and safe with high efficiency. Glutaraldehyde, in combination with O-phenylphenol and Isopropyl alcohol, provides a broad spectrum of action against pathogens such as Gram+ and Gram- bacteria, Mycobacterium tuberculosis, fungi, spores and viruses (Adenovirus, HIV, HBV and HCV). Thanks to its coformulants, it also has a strong detergent power. Barrycidal provides a rapid disinfection and allows a fast reuse of instruments. It respects the most sensitive parts (optical fibers, adhesives, plastics, rubbers and metals) and does not have a corrosive and oxidizing action.

For professional use.

Medical Device class IIb in compliance with Directive 93/42/EEC.
Instructions: GB, IT.

Made in Italy.

Technical Specifications

CHEMICAL COMPOSITION

100 g of the product contain: Glutaraldehyde 2,0 g; O-phenylphenol 0,5 g; Isopropyl alcohol; coformulants; depurated water q.s. to 100,0 g.