



Product Code	28402
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1456069
CND	Z120601
EAN/UPC	8019781501097
UMDNS	13762
PARAF	936980604

Description

T-ONE MEDISPORT ELECTROTHERAPY - 4 channels

Professional equipment specifically designed for sport specialist users. Muscle stimulation specialized programs improve sport performance. It increases strength and helps recovery of sport injuries with specific treatment protocols.

Main features

- constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- four independent channels unit (up to 8 connectable electrodes)
- compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect;
- monophasic wave generator (only for iontophoresis programs)
- 2+2 users programs
- backlight graphic colour display
- user feedback: the device recognizes patient connection
- soft-touch easy to use keyboard
- on-line 5 languages guide (GB, FR, IT, ES, DE) on display
- 4,8 V rechargeable battery pack
- personalized programs for man and woman
- maximum intensity 240 mApp, with 1 mA step regulation

Size: 150x90x35 mm

Technical Specifications

- BEAUTY (18 PROGRAMMES)
 - Firming up
 - Toning up
 - Definition
 - Modelling
 - Microlifting
 - Lipolysis abdomen
 - Lipolysis thighs
 - Lipolysis glutei and hips
 - Lipolysis arms
 - Tissue elasticity
 - Capillarization
 - Heaviness in legs
 - Face Microlifting Handy Trode
 - Breast Microlifting Handy Trode
 - Décolleté Microlifting Handy Trode
 - Face wrinkles Handy Trode
 - Forehead Wrinkles Handy Trode
 - Collagen

- REHA (10 PROGRAMMES)

Ionophoresis L
Ionophoresis M
Ionophoresis H
MENS Microcurrent
Haematomas
Oedema
Tens sequential
Tens Burst
Atrophy prevention
Atrophy

- TENS (16 PROGRAMMES)

Tens eapid
Tens Endorphinic
Tens Max Values
Inflammatory
Neck pain
Backache/Sciatic pain
Sprains/Bruises
Vascularization
Muscle Relaxant
Hand and wrist pain
Plantar stimulation
Epicondylitis
Epitroclea
Periarthritis
Tens Burst
Microcurrent

- NEMS (14 PROGRAMMES)

Warming up
Resistance
Resistance strength
Basic strength
Fast strength
Explosive strength
Deep capillarization
Muscle recovery
Agonist/Antagonist
Sequential tonic contractions
Sequential phasis contractions
Muscle Relaxant
Deep massage
EMS
30' Resistance
30' Resistance strength
20' Base strength
20' Fast strength
20' Explosive strength
20' Hypertrophy

- NEM

Tens mem 1
Tens mem 2
Tens mem 3
Tens mem 4
Tens mem 5
NEMS mem 1
NEMS mem 2
NEMS mem 3
NEMS mem 4
NEMS mem 5
NEMS alternated
Battery test

- 2+2

Tens rapid
Tens endorphinic
firming up
Toning up
Definition
Modelling
Microlifting
Muscle Relaxant
Massage
Capillarization

Inflammatory
Resistance, resistance strength
Basic strength, fast strength
Explosive strength, hypertrophy
Face Microlifting Handy Trode
Breast microlifting Handy Trode
Décolleté microlifting Handy Trode
Face wrinkles Handy Trode
Forehead wrinkles Handy Trode

Standard accessories

4 connection cables (8 connectable electrodes)
Kit of pre-gelled self-adhesive electrodes (square, rectangular, round)
Ionophoresis kit: elastic belt, 2 silicone electrodes, 2 sponges, carrying bag
Internal Ni-Mh rechargeable battery pack and battery charger
Electrode position manual and user manual (GB, FR, IT, ES). On request DE, PT