



## Disposable Biopsy Punches

- Seamless stainless steel tube blades with PP handle with and without plunger.
- Ribbed handle for sure grip and control
- Size identification on each punch
- Used for cutting, dissecting and/or filling of tissue mostly.
- The punches are sterile and single use.
- Clinical areas: Dermatology, Gynecology and Podiatry.

Sizes:


BP-10F	1.0mm	◦	BP-35F	3.5mm	○
BP-15F	1.5mm	◦	BP-40F	4.0mm	○
BP-20F	2.0mm	○	BP-50F	5.0mm	○
BP-25F	2.5mm	○	BP-60F	6.0mm	○
BP-30F	3.0mm	○	BP-80F	8.0mm	○

Length of blade: 7mm

Packaging: Each punch in a blister pack and 20 pcs. in a dispenser box, 50 boxes per export carton.

## Safe handling

- Store any new products in a dry, clean and safe place.
- Avoid high temperature, high humidity, direct sunlight, soaking and preserve at room temperature
- Do not preserve at the chemical products' place and gas occurrence place
- Preserve at the status of non-break to package material and non-pinhole
- It is necessary to confirm the use period on each package and box (Shelf life 5 years).
- By confirming the use period on each package and box, dispose the products which expire the use period.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Avoid contact of the delicate blades at tips with hard surfaces or other instruments
- Re-use prohibited
- Re-sterile prohibited
- Do not use after expiry date
- Do not use if it is considered that the products aseptic condition is spoiled by broken package.



kai industries co.ltd, 1110 Oyana, Seki City, Gifu Pref., 501-3992, Japan

kai Europe GmbH, Kottendorfer Straße 5, 42697 Solingen, Germany

CE0197

Product complies with  
Requirements of Directive  
93/42EEC