

# EASY LIFE: AIDS FOR DAILY LIVING



• **28920 CUTLERY SET**  
 • **28921 BENDABLE CUTLERY SET**  
 Set of 4 pieces (knife, fork, tea and table spoon). Made of stainless steel with thick, ribbed rubber handle (11 cm) for easy grip



• **28925 SPOON & FORK HOLDERS - pack of 2 pcs.**  
 These portable, stylish holders fit onto virtually any style spoon and fork providing a secure grip (spoon and fork are not included).



• **28966 DRINKING GLASS - box of 210**  
 Drinking glass made of polypropylene with a separated tap and a lip. Capacity 250 ml graduated by 25 ml. Autoclavable. Dishwasher safe. Not for microwaves.




Capacity: 295 ml – Diameter: 7.60 cm (3")  
 Height: 12.7 cm (5") – Weight: 0.27 kg

• **28967 HANDLE GRIP CUP**  
 Transparent plastic resistant mug with two large handles to allow a more secure grip. Extra-wide base for stability and to reduce tipping over. It includes a spouted lid with a little hole for straws. Suitable for patients with limited hand control and tremors. Mug only is dishwasher and microwave safe.




Diameter: 22.8 cm (9") Depth: 1.2 cm (1/2")

• **28930 INNER LIP PLATE**  
 High rim around the entire plate to keep food from sliding off and prevent food spills. Raised lip facilitates the grasping of food with fork or spoon. Made of durable transparent plastic, which is dishwasher and microwave safe.




Diameter: 23 cm Rim height: 2.5-5 cm

• **28931 HIGH RIM PLATE**  
 Polypropylene high rim plate with one side higher than the other to allow one-handed eating. Higher side facilitates the grasping of food with one hand without spilling it. Heat-resistant: 100-120°C Cooling resistant temperature: -20°C



Diameter: 22.8 cm (9")

• **28932 3-SECTION PLATE**  
 Plastic plate divided in 3 sections to keep food separate. It helps patients with physical disabilities, tremors, Parkinson...to eat independently. Made of resistant plastic.




Diameter: 22.8 cm (9")

• **28935 SUCTION PAD for plates**  
 White anti-slip suction pad with patented vacuum system to be attached on flat surface with slight pressure. It holds plate in place allowing patients with disease, which restrict hand movement, to eat and drink independently and reduce accidental spills.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.



• **28915 LONG BATH SPONGE**  
 Lightweight sponge features soft polyfoam attached to a plastic handle. Available with contoured sponge.



28951

• **28951 JAR OPENER & CLOSER**  
 This uniquely designed product concentrates hand power to remove lids and turn facets effortlessly.



• **28950 EASY OPENER**  
 All in one product easily opens bottle caps, pull tabs and even impossible box tops! Designed to provide safety for the user.




• **28906 MAGNIFYING GLASS**  
 Magnifies two times larger than original size. Large enough to do what's needed. Small, portable, fitting easily into purse or pocket.



Can hold 1 to 3 keys. Keys fold into handle when not in use to slip into pocket or purse.

Length 13 cm

• **28941 MULTIPLE KEY TURNER**  
 Designed to assist people with limited wrist strength or movement (arthritis or other limited hand mobility) to manipulate keys with ease. House keys and most car keys (without plastic head covers) fit into the slots and are held by a single screw.



With cover

• **28912 SOCK AID NOTCH**  
 Soft terry cloth side grips the sock, while smooth nylon side reduces friction, allowing foot to slip in easily. Two long loop handles are easy to hold and assist with pulling the sock onto foot. Ideal for people who can use one hand only or with decreased strength or coordination.



28911

• **28911 ZIPPER AND BUTTON PULLER**  
 It helps pulling the button or zipper into position. Folds to half-length to fit easily into pocket or purse. Ideal for people with weak or arthritic fingers.



• **28945 SHOPPING BAG CARRIER**  
 This product provides an effort-free handle that remains constantly supportive, regardless of bag loads.