

GRADUAL COMPRESSION, UNISEX, COTTON SOCKS

UNISEX GRADUAL COMPRESSION COTTON SOCKS

High quality knee socks made of soft cotton combined with elastic fibres to provide proper gradual compression. Special attention was applied to the development of this cotton product line to ensure maximum fibre quality, along with a soft elastic cuff that offers excellent levels of expandibility: all to ensure maximum comfort and wellness throughout the day.

This product line is available in two compression levels (medium and strong) and in two colours (black and blue). On request, for quantities of 100 or more they are also available in light and extra strong compression and beige and grey colours.

Medium compression level: 14-16 mm/Hg at the ankle

To help prevent discomfort and heaviness in the legs.
Composition: 67% cotton, 26% polyamide, 7% elastane.

High compression level: 15-21 mm/Hg at the ankle

To help prevent the onset of venous disease.
Composition: 63% cotton, 27% polyamide, 10% elastane.
It is a Medical Device.



MADE IN ITALY



Multilingual box and instructions:
GB, FR, IT, ES, PT, DE, GR, HU



GRADUAL COMPRESSION COTTON SOCKS - pair

MEDIUM COMPRESSION SOCKS

Colour	Small	Medium	Large	X-Large	XX-Large
Black	41220	41221	41222	41223	41224
Dark blue	41225	41226	41227	41228	41229

STRONG COMPRESSION SOCKS

Black	41240	41241	41242	41243	41244
Dark blue	41245	41246	41247	41248	41249

Shoe size

European	37 - 39	39 - 41	41 - 43	43 - 45	45 - 47
US - Men	5 - 6½	6½ - 8	8 - 9½	9½ - 11	11 - 12½
US - Women	6½ - 8	8 - 10	10 & over	-	-

SUPPORT SHEER LINE

Knee Socks, Stockings and Tights designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and super soft fabrics.

Available in two compression levels (medium and strong) and in beige or black.

Medium compression level: 70 den - 10-14 mm/Hg at the ankle

Ideal for those who would like to initiate a prevention programme, but are mainly interested in the comfort and well-being offered by elastic compression stockings.

Composition: 84% polyamide, 16% elastane (Stay ups, Pantyhoses).
85% polyamide, 15% elastane (Knee-highs).



Multilingual box and instructions:
GB, FR, IT, ES, PT, DE, GR, HU

MADE IN ITALY

Knee-high with a gradual compression,
comfort top, honeycomb mesh



Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, HU



MADE IN ITALY

Stay ups with gradual compression,
silicone lace top, honeycomb mesh



Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, HU, NL

MADE IN ITALY

Pantyhose with gradual compression,
reinforced body, honeycomb mesh

KNEE-HIGHS - pair

MEDIUM COMPRESSION

Colour	Small/Medium	Large/X-Large
Black	41100	41101
Beige	41102	41103

STRONG COMPRESSION

Black	41106	41107
Beige	41108	41109

STAY UPS - pair

MEDIUM COMPRESSION

Colour	Small	Medium	Large	X-Large
Black	41120	41121	41122	41123
Beige	41124	41125	41126	41127

STRONG COMPRESSION

Black	41132	41133	41134	41135
Beige	41136	41137	41138	41139

PANTYHOSES

MEDIUM COMPRESSION

Colour	Small	Medium	Large	X-Large	XX-Large
Black	41150	41151	41152	41153	41154
Beige	41155	41156	41157	41158	41159

STRONG COMPRESSION

Black	41165	41166	41167	41168	41169
Beige	41170	41171	41172	41173	41174