


## RAISED TOILET SEATS




• **27750 RAISED TOILET SEAT** height 12.5 cm with fixing system (front screw) and fixed armrest made of anodized aluminium. Max load 130 kg (286.6 lb). Packing: 410x165x468 mm



• **27741 RAISED TOILET SEAT** - 10 cm with lateral fixing system. Max load 250 kg (551 lb). Packing: 373x128x410 mm



• **27751 RAISED TOILET SEAT** height 12.5 cm with fixing system (back and front screw) and anodized aluminium reclining armrest. Max load 130 kg (286.6 lb). Packing: 410x165x468 mm



• **27740 RAISED TOILET SEAT** - 15 cm with lateral fixing system. Max load 250 kg (551 lb). Packing: 373x165x410 mm




• **27756 CLIPPER RAISED TOILET SEAT** - 11 cm Easily adapts to any toilet bowl. Totally made of injected plastic. 4 fixation clips to secure strongly the seat on the toilet bowl. Max load 185 kg (407 lb).



• **27735 RAISED TOILET SEAT with LID** - 10 cm  
• **27737 RAISED TOILET SEAT with LID** - 15 cm with lateral fixing system. Max load 250 kg (551 lb). Packing: 373x128x410 mm (27735) 373x165x410 mm (27737)

## CONTACT PLUS SOFT RAISED TOILET SEATS

• **27743 CONTACT PLUS SOFT RAISED TOILET SEAT**  
• **27745 CONTACT PLUS SOFT RAISED TOILET SEAT - with lid**  
High quality soft toilet raiser, made of high-density foam injected in a waterproof PVC skin. It fits any toilet bowl in a very secure way. This comfortable seat with ergonomic shape and high density foam, raises the toilet by 11 cm (4.3 inches). Tested for 185 kg (407 lb).



Dimension	width	depth	height
Outside	37 cm	41 cm	16 cm
Inside	21 cm	26 cm	16 cm




## WC - SHOWER COMMUNE



• **27758 BATH BENCH with backrest** Aluminium frame and padded seat. Removable plastic bowl WC.

• **27759 WC - SHOWER COMMUNE** Suitable also as raised toilet seat. Adjustable height 45-55 cm (back 78-88 cm)  
• **27766 PLASTIC WC BOWL - spare**

## CLIPUP® TOILET RAISER



MADE IN FRANCE



• **27738 CLIPUP® TOILET RAISER** - autoclavable  
ClipUp® can be installed safely on any type of toilet bowl. Innovating fixation system with ratchet. Design with front cutout provides optimal comfort and facilitates personal hygiene. Max Load 150 kg.

TECHNICAL SPECIFICATIONS			
Size:	width 34 to 37 cm	depth 38.5 cm	height 14 cm
Packing:	45.5x16.5x 39.5 cm	Seat opening (depth x width): 26x16 to 19 cm	Toilet elevation: 10 cm
		Weight: 0.872 kg	Material: 100% polypropylene