

## COMMODORE WHEELCHAIRS

Packing size:  
54.5x19x76 cm  
Weight: 7.7 kg



43196

4 castors Ø 10 cm  
(4 ")

- 43196 BLUE COMMODORE WHEELCHAIR
- 43198 SMART COMMODORE WHEELCHAIR

Made of aluminium with anodized surface. Fixed height. Equipped with 4 castors with brakes, one large footrest, and pan with lid. Indoor use. Multilingual manual: GB, FR, IT, ES, PT, DE.

Packing size:  
72x58x36 cm  
Weight: 8.8 kg



43198

4 castors Ø 10 cm  
(4 ")

GIMA code						
43196	82	92	56.5	51	46	45
43198	96	100	60	53	52	50

## ELITE COMMODORE



43203

Size:  
54x88xh 105 cm  
Weight: 21 kg  
Load: 100 kg

43203

Supplied with plastic pan



4 wheels Ø 20 cm (2 with brake)

- 43203 ELITE 200 COMMODORE - 43 cm
- 43204 ELITE 200 COMMODORE - 53 cm
- chromed steel frame
- blue, washable, flame retardant leatherette
- adjustable length and raising footrest
- removable armrest - reclining backrest

## COMMODORE CHAIRS

- 27701 COMMODORE CHAIR - chromed steel
- Chromed steel chair frame, with plastic detachable commode seat.

- 27702 COMMODORE WHEELCHAIR - painted
- Powder coated steel frame, detachable armrest and footrest, solid rear castors with lock, plastic commode seat.

- 27703 COMMODORE WHEELCHAIR - stainless steel

Stainless steel chair frame, detachable armrest and footrest, rear castors with lock, plastic commode seat. Multilingual manual: GB, FR, IT, ES.

- 27699 FOOTREST

for 27702 - pair

- 27700 FOOTREST

for 27703 - pair

- 27766 PLASTIC WC BOWL - spare - Ø 26xh 15 cm



27701



27702



27703

GIMA code												
27702/3	55	44	12.5	12.5	51	88	90	44	42	100	54x53x53	17 18
27701	56	44			47	94		43	38	100	50x19x75	8.5

## COMMODORE AND SHOWER CHAIR

- 43193 COMMODORE AND SHOWER CHAIR

3 in 1 aluminium lifting chair which can be used as raised toilet seat, commode chair and shower chair. Designed with spring-loaded uplift mechanism and safety lock, that help users to stand up, pushing on the armrests. Self-powered, lean forward armrests for easy exit. Anti-slip tips with height adjustable legs.

Delivered with bucket, bucket cover and splash-shield. Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.

### TECHNICAL SPECIFICATIONS

Width: 52-55.8 cm  
Depth: 69.2-73 cm  
Height: 85-95.5 cm  
Packing: 64.5x58.5x39.5 cm  
Max load: 136 kg - Weight: 7.1 kg

Cushion and backrest made of waterproof durable quality PU foam

3 IN 1:  
COMMODORE CHAIR,  
SHOWER CHAIR  
TOILET RAISER



43193

Adjustable seat height

Spring-loaded uplift mechanism

Self powered, lean forward armrest for easy exit

Rugged, rust-proof aluminium frame

## COMFORT COMMODORE

- 43190 COMFORT COMMODORE CHAIR - blue

- 43191 COMFORT COMMODORE CHAIR - height adjustable - blue

Padded and removable seat for best comfort and high level of hygiene. Pan is removable from both sides and from the top.

MADE IN FRANCE

MAX LOAD 130 kg

### TECHNICAL SPECIFICATIONS

Total height: 85.5 or 85.5/97 cm  
Width between armrests: 56 cm  
Total width: 63 cm  
Height armrest-floor: 71 or 71/83 cm  
Height seat-floor: 50 or 50/63 cm  
Height seat-backrest: 40 cm  
Seat plate depth: 39 cm  
Total depth: 49 cm



43190

Removable seatplate

Weight: 6.6 kg

Comfortable padding