

SKID-PROOF STICKS WITH LIGHT

- 43044 TRUSTY CANE - silver
- 43045 TRUSTY CANE - black

Sturdy, folding cane with built-in lights and extra wide pivoting head with 3 firm and skid-proof bases. Height adjustable 85-97 cm. Large base with Ø 9 cm. Shaft in aluminium (thickness 1.2 mm) Handle in ABS (15 cm). Weight 420 g. Need 3x1.5 V AAA batteries not included.

- 43048 FOOT PADS - kit of 3 - spare



STANDS ON ITS OWN

Extra wide pivoting head revolves 360° degree



STICKS WITH LIGHT AND ALARM



- 43050 SAFETY LIGHT STICK - black
 - 43052 SAFETY LIGHT STICK - skid proof base - black
- Innovative sticks with safety light and alarm.
- front ultra-bright white 4 LED light up to 15 lumen, with adjustable positions
 - ultra-bright red LED warning flash with 2 LED, up to 5 lumen
 - ABS hand grip and aluminium frame
 - 100 dB warning siren
 - adjustable height: 10 positions from 78 to 100.5 cm
 - extra wide 360° pivoting base (43052 only)
 - stands on its own (43052 only)
- Works with 2 AAA batteries, not included. Maximum load: 110 kg

Supplied with carrying hand strap

43050 43052

QUADRIPODS AND STICK WITH SEAT



- 43060 STICK WITH SEAT - silver
- Aluminium walking stick with seat, height 84 cm, weight 0.9 kg
- 43065 TRIPOD - silver
- Adjustable height 74-97 cm, weight 0.8 kg
- 43070 QUADRIPOD Q1 with PU handle - aluminium - silver
- Adjustable height 74-97 cm, weight 0.9 kg
- 43071 QUADRIPOD Q2 with small base and PU handle - silver
- Adjustable height: 65-88 cm, weight 0.8 kg
- 43072 QUADRIPOD Q3 with PU handle - bronze
- Adjustable height: 70-88 cm - weight 0.5 kg
- Maximum load for tripod and quadripods: 100 kg
- 27796 RUBBER FERRULE int. Ø 16 mm - grey for 43065, 43070 - bag of 5



CANE, TRIPOD, QUAD CANE



- CANE, TRIPOD AND QUAD CANE WITH ADJUSTABLE HEIGHT
- Walking aids for people who need to relieve the weight from the legs providing stability, higher than a traditional stick.
- 43075 STABILE WALKING CANE
 - 43076 TRIPOD WITH RING HANDGRIP
 - 43077 QUAD CANE WITH RING HANDGRIP
- Stick Structure: all the bases and the upper sliding tube are in light aluminium alloy. For the quad cane the central tube is made of steel. U-shaped tubes with strong injection-moulded plastic joint to provide strength and stiffness of the base (only 43076-7).
- Tube Locking clip: locking of the telescopic tubes with a clip equipped with a double pin; this improves the safety.
- Tapered ring nut: supplied with a tapered ring nut allowing, by tightening, to eliminate noise.
- Tips: high grip and have Ø of 35 mm, ISO EN 11334-4. Tip made of high friction material (only 43075)
- Class I Medical Devices 93/42/EEC.