

WALKING FRAMES



• 27721 WALKING FRAME - 2 WHEELS

Intended to provide support to the patient. It allows mobility making it easy to walk.

Size 54x57 cm with adjustable height 52-87 cm.

Instructions: GB, FR, IT, SE.

Weight: 8 kg - Load: 100 kg.

• 27722 WALKING FRAME

Made of lightweight aluminium tubing, non-slip end and padded handgrips.

Size 47x60 cm foldable with adjustable height 78-96 cm. Instructions: GB, FR, IT, SE.

Weight: 2.4 kg - Load: 100 kg.

• 27794 RUBBER FERRULES Ø 28 mm - grey for 27722-3 - bag of 5 pcs.



WALKING FRAMES WITH WHEELS




• 27723 WALKING FRAME - 2 WHEELS

Made of lightweight aluminium tubing, with padded hand grips, 2 non-slip ends and 2 with wheels. Size 51x64 cm, foldable with adjustable height 79-97 cm.

Instructions: GB, FR, IT, SE.

Weight: 3 kg - Load: 100 kg.

• 27794 RUBBER FERRULES Ø 28 mm - grey for 27722-3 - bag of 5 pcs.



• 43129 WALKING FRAME - 2 INTERNAL WHEELS


Made of lightweight aluminium tubing with padded handgrips, 2 non-slip ends and 2 with wheels. Foldable.

Adjustable height 80-90 cm.

Instructions: GB, FR, IT, SE.

Weight: 2.5 kg - Load: 100 kg.

Size 51x58 cm.



• 43130 WALKING FRAME - 4 WHEELS

Made of lightweight aluminium tubing, 2 front swivel wheels, 2 rear wheels with stop, with foam grips.

Size 55x65 cm foldable with adjustable height 83-93 cm.

Instructions: GB, FR, IT, SE.

Weight: 3.5 kg. Load: 100 kg.

AXILLARY WALKING FRAME



• 43134 AXILLARY WALKING FRAME - 4 WHEELS - blue

Walking frame made of steel with 4 wheels and 2 brakes. Equipped with removable seat pad and soft PU foam arm pad.

Height adjustable armrests 114-139 cm.

Height adjustable back support 95-110 cm.

Foldable: 15.5 cm width.

Frame size: 71x54.5x h 71 cm.

Max load: 120 kg.

Weight: 10.5 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.

STAND UP ROLLATOR

• 43172 FULLY ADJUSTABLE STAND UP ROLLATOR

It helps adults to stand upright and look ahead with better support.

Compact one hand folding design, aluminum stable frame, dual hand brake system and fully adjustable armrest, backrest and length.

Detachable front wheels and stand up handles.

Padded seat and backrest for more comfort.

Two front (Ø 25.4 cm) and two rear (Ø 20.3 cm) PVC wheels,

cane holder and shopping bag.

Packing size: 78x28.5xh 45 cm

Weight: 12 kg. Load: 136 kg.

