

☐ ELECTRIC WHEELCHAIRS - FOLDABLE



- **43290 ELECTRIC WHEELCHAIR** - 30 cm rear castors
- **43291 ELECTRIC WHEELCHAIR** - 56 cm rear wheels with handrail

Electric wheelchair, for daily indoor and outdoor use. It can be used in electric or manual mode, for easy clutch switching. User can move wheelchair forwards, backwards or turn around through an intelligent controller, complete with joystick to control speed and direction. It is provided with aluminium alloy frame, hand rest, foot rest, expandable arm rest and flexible shock-proof suspension. Removable seat and backrest cushions, battery and remote control.

Comfortable driving with ergonomic design which relieves fatigue.

Powerful: efficient energy-saving motor; high-precision gearbox; fast and stable starting. 5 speed settings.

Triple protection: electromagnetic brake and gear self-lock protection.

Stylish model: smooth light design characterized by dynamic guard plate, split fixed seat and patented arc-shaped foot rest. Backrest folding mechanism.

Supplied with lithium battery and charger, manual (GB, IT).

- **43292 LITHIUM BATTERY** - spare for 43290-1
- **43293 PNEUMATIC FRONT CASTOR 20 cm (8")** - spare for 43290-1
- **43294 PNEUMATIC REAR CASTOR 30 cm (12")** - spare for 43290
- **43295 PNEUMATIC REAR WHEEL 56 cm (22")** - spare for 43291




TECHNICAL SPECIFICATIONS

Load capacity:	≤100 kg
Enduranced mileage:	≥20 km
Maximum speed:	≤6 km/h
Battery:	lithium battery, DC 24V x 18Ah
Motor power:	130 W
Power consumption:	≤2.5 kw x h/100 km
Horizontal brake performance:	≤1.5 m
Maximum safe ramp brake:	≤3.6 m (6°)
Stability:	≥9° (static), ≥6° (dynamic)
Noise:	< 65 dB
Obstacle crossing height:	≥40 mm
Minimum turning radius:	1.2 m
Groove clearance width:	100 mm
Hill-holding performance:	9°
Size:	670 mm - open, 370 mm - closed (43290) 660 mm - open, 360 mm - closed (43291)
Seat width:	440 mm
Net weight:	26 kg (43290) - 29 kg (43291)
Protection:	IPX4


☐ PEDAL EXERCISERS

Compact exercise bikes for moderate cardiovascular workout. The regular use stimulates blood circulation and increases

muscle strength. Can be placed on a tabletop to exercise arms or on the floor as a bicycle-type exerciser while sitting.



- **43150 BASIC PEDAL EXERCISER**
Delivered in kit with screw and key for easy assembling.



- **43151 DISPLAY PEDAL EXERCISER**
Foldable to save space. 5 function display: exercise time, rotation count, rotations per minute, calories burned and scan.



- **43153 ELECTRIC PEDAL EXERCISER WITH DISPLAY**
Made of ABS and steel with integrated LED display to track and read time, speed and calories. 6 adjustable speed levels. Infrared remote control (battery not included). Power supply: 80 W, 220V-50 Hz.