

UNIVERSAL BELTS FOR ARMCHAIRS, COMMODOES, WHEELCHAIRS

SMART LINE (43216, 43217)

MADE IN ITALY

Padded harness and belt to prevent patient sitting on wheelchair from sliding forward. Equipped with two black propylene webbing (height 5 cm), adjustable up to 100 cm. One-hand release safety clip closure. Harness has two adjustable suspenders with sturdy strap. Class I, Medical Device. Washable at 30°. Harness size: 60xh 37 cm, max. adjustable length: 260 cm. Belt size: 48xh 15 cm, max. adjustable length: 248 cm.



• 43216 ABDOMINAL BELT



• 43217 HARNESS

PREMIUM LINE (43175, 43176, 43177)

MADE IN ITALY

Aids used for positioning and handling patients with safety. Equipped with button acetate closure to increase safety level and prevent it from releasing. Ideal for wheelchairs, commodes and armchairs (with 43178).

Coating in class 1 fire-retardant and water-repellent fabric. Washable at max. 40°

• 43178 EXTENDER FOR PREMIUM

LINE BELTS - up to 150 cm

Adjustable length.

Necessary for use with armchair.



• 43175 PELVIC BELT



• 43176 ABDOMINAL BELT



• 43177 HARNESS



BED RESTRAINT BELTS

100% cotton safety belts used for control of bedridden patient. Equipped with 1 (43185) or 3 (43186-7) safety locks that can be released with a special magnetic key. The belts allow a single (43185-7) or double (43186-7) movement to position the patient in a supine (43185-7) or lateral position (only 43186-7). Washable at 90° (43185) or 95° (43186-7). Belt height 21 cm.

Safety locks and magnetic key



43185

MADE IN ITALY

GIMA code	BED RESTRAINT BELTS	Size	Circumference cm
43185	Single movement belt	L-XL	80-122
43186	Double movement belt	S-M	60-88
43187	Double movement belt	L-XL	80-122

43186, 43187



TRANSFER BELT

• 43181 TRANSAC TRANSFER BELT - small

• 43182 TRANSAC TRANSFER BELT - medium

Polyurethane transfer belt to be placed around the patient's or the caregiver's waist, to help the patient sit or stand with minimum effort and maximum safety. It has 5 vertical and 4 horizontal padded handles with polyester strap and nylon loop. Length: 105 cm (S), 120 cm (M). Transferring load: 160 kg. Washable at 60°.

43182



MADE IN FRANCE

ROTATING SEAT CUSHION

• 43040 ROTATING SEAT CUSHION

Cushion to assist someone to turn easily in a car, on the bed or on the chair. Size: Ø 40 cm

43040



Slippery nylon filled with PU foam

Upper in 100% polyester




Bottom in non-slip PVC fabric

PVC WATER MATTRESS

• 28502 PVC WATER MATTRESS

It helps movement and distributes pressure on the whole supporting surface. Double sealed seams provide maximum strength and durability. Highly resistant PVC material construction easy to clean with any mild detergent. Air frame + 3 water elements. Supplied with hand pump and hose.



28502

TECHNICAL SPECIFICATIONS

Material: PVC 0.40 mm thickness
Size of outer air frame: 200x93x11.5 cm
Size of each element: 69x60x11.5 cm
Weight: 3.5 kg

SHAMPOO BASIN

• 28507 SHAMPOO BASIN

Inflatable shampoo basin ideal for personal in-bed hair care needs. It supports the head, neck and shoulders during shampooing. A 20 cm deep basin helps prevent splashing. Drain hose efficiently removes water and foam. Supplied with water bag with shower head in clip. Size: 53x61x20 cm Capacity: 10 l

• 28509 WATER BAG - spare

28509

Bag for easy water transport when water source is not easily accessible



28507