

LIFT ARMCHAIRS WITH 1 OR 2 MOTORS

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position. Armchairs and covers made of microfibre, easy to clean. Power supply: 100-240 V - 50/60 Hz. Multilingual manual: GB, FR, IT, ES.

1 motor lift armchair:

the chair's backrest and footrest can lift or recline together.

2 motors lift armchair with extra lifting footrest:

it can reach the bed position providing a relaxing experience. Its backrest and footrest can operate separately to any position. Extra lifting footrest allows 8 cm lifting adjustment higher than flat position.

ARIANNA LIFT ARMCHAIRS with 1 or 2 motors

Blue	Brown	Green	Burgundy	Motors	Fireproof	Footrest
45170	-	-	45173	1 motor	no	standard
45175	45176	45177	45178	2 motors	no	extra
45179	-	-	-	2 motors	yes	extra

COVERS FOR LIFT ARMCHAIRS

45195	45196	45197	45198	Covers for armchairs
-------	-------	-------	-------	----------------------

SPARE MOTORS FOR LIFT ARMCHAIRS

45162	Motor 1 for 1 motor chair - spare
45163	Motor 2A for 2 motors chair - up and down - spare
45164	Motor 2B for 2 motors chair - back and footrest - spare

Seat			Backrest	Chair		Load	
Width cm	Depth cm	Height cm	Height cm	Width cm	Length cm (reclined)	kg	lbs
49	49.5	49	70	75	170.5	150	331

ARIANNA LIFT ARMCHAIRS 1 or 2 MOTORS



45176



45178

45175,
45179

45183

GINEVRA LIFT ARMCHAIRS 2 MOTORS



45180



45182

GINEVRA LIFT ARMCHAIRS with 2 motors

Blue	Brown	Green	Burgundy	Motors	Fireproof	Footrest
45180	-	45182	45183	2 motors	no	extra
-	45184	-	-	2 motors	yes	extra

COVERS FOR LIFT ARMCHAIRS

45195	45196	45197	45198	Covers for armchairs
-------	-------	-------	-------	----------------------

SPARE MOTORS FOR LIFT ARMCHAIRS

45163	Motor 2A - up and down - spare
45164	Motor 2B - back and footrest - spare

Seat			Backrest	Chair		Load	
Width cm	Depth cm	Height cm	Height cm	Width cm	Length cm (reclined)	kg	lbs
46	51	40	66	74	169	150	331

FLAVIA LIFT ARMCHAIRS with 2 motors

Blue	Brown	Burgundy	Motors	Fireproof	Footrest
45185	45186	45188	2 motors	no	extra
-	-	45190	2 motors	yes	extra

COVERS FOR LIFT ARMCHAIRS

45195	45196	45198	Covers for armchairs
-------	-------	-------	----------------------

SPARE MOTORS FOR LIFT ARMCHAIRS

45163	Motor 2A - up and down - spare
45164	Motor 2B - back and footrest - spare

Seat			Backrest	Chair		Load	
Width cm	Depth cm	Height cm	Height cm	Width cm	Length cm (reclined)	kg	lbs
47	50	47	74	78	171.5	150	331

8 CM EXTRA LIFTING FOOTREST

Extra lifting footrest allows 8 cm lifting adjustment higher than flat position



FLAVIA LIFT ARMCHAIRS 2 MOTORS

45188,
45190

45185



45186