

## HOSPITAL BEDS



## BED WITH ADJUSTABLE HEAD SECTION

## • 27600 BED WITH 1 ARTICULATION

## • 27601 BED WITH 1 ARTICULATION with castors

Painted steel frame with net.

Pipe frame with net for head movement.

2 heads in chromed pipe complete with MDF bilaminated coated panel. Movement by a crank. Protection of end run.

Rubber castors Ø 100 mm, with brakes (only for code 27601)

Size: 206x89xh 90 cm. Height of frame from bottom 52/55 cm.

Size of frame: 200x85 cm. Weight 32 kg. Load 110 kg.



## • 27607 STANDARD PLUS BED

## • 27608 STANDARD PLUS BED with Ø 100 mm castors

Chrome-plated steel bed ends.

Laminated plastic panels. Rigid mesh mattress platform.

Adjustable back-rest. Size: 206x89xh 90 cm. Load 110 kg.

Height adjustable  
160 to 245 cm

MADE IN ITALY



## BED WITH ADJUSTABLE HEAD AND FOOT SECTION

## • 27610 BED WITH 3 ARTICULATIONS

## • 27611 BED WITH 3 ARTICULATIONS with castors

Painted steel frame with 3 section net.

Painted steel pipe, 200x85 cm.

2 heads in chromed pipe complete with MDF bilaminated coated panel. Movement by 2 cranks. Protection of end run.

Rubber wheels Ø 100 mm, two with brake (only 27611)

Size 206x89xh 90 cm. Height of frame from bottom 52/55 cm.

Size of frame: 200x85 cm. Load 110 kg. Weight 42 kg.



## • 27605 STANDARD BED

## • 27606 STANDARD BED with Ø 50 mm castors

Basic bed made of epoxy - coated steel.

Head and foot ends with plastic feet (27605) or castors (27606).

Rigid mesh mattress platform. Load 110 kg.

## ACCESSORIES FOR BEDS



## • 27697 N-SHAPED BEDRAIL 165x40 cm

Bedrail made of steel frame and surface with coating.

Siderail can be fixed by screw. Max load: 75 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, RO, SE, GR.



## • 27698 DELUXE ADJUSTABLE BEDRAIL 194-220x43.5 cm

Bedrail made of steel frame and surface with coating.

Adjustable siderail length by spring. Max load: 75 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, RO, SE, GR.



## Bed side pad

Protection for patients against accidental collisions on rigid sides of the bed.

Made of polyurethane foam with not removable outer coating water-repellent and flame retardant.


Size: 166x50xh 3 cm



## Raising block

Designed to raise bed or chair by 9 cm (inner

Ø 8.7 cm, base Ø 16.6 cm)



## IV stand

Stainless steel IV stand with 4 hooks

GIMA code	ACCESSORIES FOR ALL BEDS
27641	Universal collapsible supports - pair - see above
27644	Lifting pole - see above
27646	4 Hook IV stand for bed - stainless steel
27637	Bed side pad with velcro
27649	Raising blocks 9 cm - kit of 4
27841	I.V. Stand - 5 wheels trolley, see page 140