

■ DRUG TEST READER - FOR URINE DRUG CUPS

- 24561 DRUG CUP READER with software and printer
- 24562 MULTIDRUG CUP - 7 drugs
- 24563 MULTIDRUG CUP - 8 drugs

Gima professional drug test reader digitally scans DOA urine cups and gets a complete analysis on touch screen.

It displays and prints the qualitative results only.

- automatic test menu recognition (2D barcode)
- in-process system quality control check
- memory card to store test data
- data management/export, archived results

Main benefits:

- automation: it eliminates results interpretation variability/subjectivity
- accurate: performance advantages vs. visual interpretation
- user-friendly: fast results and automatic print
- portable design
- LIS connectivity

Built-in printer and software in GB, FR, IT, ES, PT, DE, GR, HR.

Supplied with print paper roll, power cord and adapter, EU adapter, calibration cup, touch pen and user manual in GB, FR, IT, ES, PT, DE, GR, HR.

- 35674 THERMAL PAPER - 57 mm x 25 m - spare

24562

USB and
LAN portTouch pen
included

TECHNICAL SPECIFICATIONS

Mode:	photoelectric conversion
Display:	7" LCD touch screen
Reading time:	1 - 2.5 minutes
Detection type:	qualitative
Light source:	LED
Connection:	LAN port, 4 USB ports
Printer:	built-in thermal printer
Power supply:	DC 12V, 2A
Operating altitude:	0-2000 m
Size:	221x215xh 165 mm
Weight:	< 1.5 kg

DRUG	24562	24563	CUT OFF NG/ML	Detection time
MTD	•	•	300	5-7 days
MDMA	•	•	500	5-7 days
MET	•	•	500	3-5 days
COC	•	•	300	2-5 days
THC	•	•	50	-
AMP	•	•	500	2-5 days
MOP	•	•	300	5-7 days
BUP	•	•	10	5-7 days
S.V.T.	•	•		
OX	•	•		
S.G.	•	•		
PH	•	•		
NiT	•	•		
Glut	•	•		
CRE	•	•		

24561

■ MULTI DRUG URINE TEST CUP - PROFESSIONAL

- 24545 MULTI DRUG URINE TEST CUP - 7 drugs - box of 25 tests

A rapid test for simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.

Detects following parameters: (MTD/MDMA)+<COC>+<AMP500>+<MET500>+<THC>+<MOP>.

- easy to use: no urine handling, just insert the key into the active test and read results in 5 min.

- tightly sealed urine cup: no leaking or spilling of urine specimen

- special design: ensures that urine used in the test, is separated from residual urine, that will remain uncontaminated, to be used for a further test

- key separately packed: only technician can take the test.

Professional use. Manual: GB, FR, IT, ES, PT, DE.

Precision: the results gained $\geq 75\%$ accuracy in $\pm 25\%$ cut-off level specimen and 100% accuracy in negative

SCREENING TEST CUP - 24545

DRUG	CUT OFF NG/ML	DETECTION TIME*
MOP	300	5-7 days
COC	300	2-5 days
THC	50	-
AMP	500	2-5 days
MDMA	500	5-7 days
MTD	300	5-7 days
MET	500	3-5 days

*Drug use over detection time not measured



Test panel
is tightly
fastened
to the cup,
protecting
it from
contamination

Max shelf life
22 months

■ MULTI-DRUG (6 DRUGS - 8 PARAMETERS) MIDSTREAM SALIVA RAPID TEST

- 24566 MULTI DRUG SALIVA MIDSTREAM TEST - box of 20 tests

A multi parameters, all-in-one test for simultaneous detection of 6 different drugs in 10 minutes. Cost effective system for the screening of workers for illicit substance consumption.

Professional use. Multilingual box and manual: GB, FR, IT, ES.

24566



Drug test	DRUG OF ABUSE SALIVA MULTITEST	Parameters	Cut-off (ng/ml)
Opiates (OPI/MOP):	Morphine		40
	6-Monoacetylmorphine		4
Cocaine (COC):	Benzoylcegonine		30
Amphetamine (AMP):	d-Amphetamine		40
Oxycodone (OXY):	Oxycodone		40
Methamphetamine (MET):	d-Methamphetamine		40
	3, 4-Methylenedioxyamphetamine		50
	(MDMA)		
Marijuana (THC):	11-nor- Δ^9 -THC-9 COOH		10


■ DRUGS OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

- 24540 COCAINE TEST - strip - urine - cut-off: 300 ng/ml
- 24541 MARIJUANA TEST - strip - urine - cut-off: 50 ng/ml
- 24542 OPIATES TEST - strip - urine - cut-off: 2.000 ng/ml

Codes 24540-41-42 are sold in box of 50 strips individually packed in single pouch.



24540-41-42



24550

- 24550 MULTI DRUG PANEL TESTS - urine - box of 25 tests for 6 drugs*
- 24549 MULTI DRUG PANEL TESTS - urine - box of 25 tests for 8 drugs**
- 24551 MULTI DRUG PANEL TESTS - urine - box of 25 tests for 12 drugs*

Accurate results in 5 min. on the following drugs: amphetamine*, cocaine*, marijuana cannabis*, methamphetamine*, opiates*, phencyclidine*, barbiturates, benzodiazepines, ecstasy, morphine300, methadone, tricyclic antidepressants.

*24550 contains tests for drugs with asterisk

**24549 contains tests for drugs with asterisk plus MDMA and 6-AM.

Cut-off levels meet SAMHSA guidelines.
Manual: GB, FR, IT, ES, DE, PT, GR, Arabic.

24551



Drug	24550	24549	24551	Detection time
MOP	-	-	300	5-7 days
COC	150	150	150	2-5 days
THC	50	50	50	5-8 days
AMP	500	500	500	2-5 days
MDMA	-	500	500	5-7 days
MTD	-	-	300	5-7 days
MET	500	500	500	3-5 days
BZO	-	-	300	7-10 days
BAR	-	-	300	2 days - 4 weeks
PCP	25	25	25	-
OPI	2000	2000	2000	5-7 days
TCA	-	-	1000	-
6-AM	-	10	-	-