

ALUMINIUM TROLLEYS



- **27429 GIMA 2 TROLLEY**
- small: 60x40xh 77 cm - weight 11 kg
 - **27425 GIMA 2 TROLLEY**
- medium 70x50xh 78 cm - weight 14 kg
- Practical trolleys with aluminium body and handles. Two white bilaminated shelves 18 mm thick, black PVC bordered, with guard-rail. Two large lateral handles facilitate transport or lift up of the trolley. 4 revolving nylon castors Ø 80 mm. Max load (x shelf): 40 kg.
- **27426 GIMA 3 TROLLEY** - medium 70x50xh 78 - weight 15 kg. Same as GIMA 2, but without guard rail. Delivered in kit form.

ANODIZED ALLOY TROLLEYS



- **27460 DELUXE TROLLEY**
 - **27461 DELUXE TROLLEY with drawer**
 - **27422 DELUXE TROLLEY with 2 drawers**
- Size: 65x45xh 80 cm
Drawer size: ext. 58x40xh 10.5 cm, int. 51.5x38.5xh 8.5 cm.
Top quality trolleys made of lightweight anodized special alloy with detachable stainless steel shelves.
4 revolving rubber castors Ø 80 mm.
Max load (x shelf): 20 kg.
Instructions: GB, FR, IT, ES PT, DE, GR.
Delivered in kit form.

LANCART: STAINLESS STEEL TROLLEYS TO BUILD UP



LANCART TROLLEYS

Stainless steel trolleys mounted on castors Ø 80 mm. Maximum trolley load: 60 kg. Range of optional accessories to customize your trolley. Drawer size (45808-45809 only): 29x30xh 7.5 cm

MADE IN ITALY

Delivered in kit form.

GIMA code	LANCART TROLLEYS	Size	Shelves	Drawer
45810	Lancart trolley - small	60x40xh 80	2	-
45811	Lancart trolley - medium	72x50xh 80	2	-
45812	Lancart Ter trolley - small	60x40xh 80	3	-
45813	Lancart Ter trolley - medium	72x50xh 80	3	-
45808	Lancart trolley - small	60x40xh 80	2	yes
45809	Lancart trolley - medium	72x50xh 80	2	yes
GIMA code	ACCESSORIES FOR LANCART			
45815	Set of siderails (for 1 shelf) for 45808, 45810, 45812			
45816	Set of siderails (for 1 shelf) for 45809, 45811, 45813			
45817	Bottle rack			
45818	Waste bin			

