


SMART CARTS - WITH 2, 3 OR 4 SHELVES OR DRAWER



SMART CARTS - grey/white RAL 7035
Carts with base and 2, 3, or 4 wide, sturdy shelves in 3 different widths 40, 50 or 60 cm.
Aluminium columns and 4 twinned castors Ø 65 mm (2 with brake).
First shelf can be mounted flat or tilted $\pm 10^\circ$.
Fully removable drawer (27893-4) without ball bearing guides.
Designed for medical devices.
Illustrated assembly instructions.

TECHNICAL SPECIFICATIONS

Dimensions (cm): 48,2x50x78 (27893, 27895)
58,2x50x78 (27894, 27889)
48,2x50x88 (27896)
48,2x50x115 (27897)
68,2x50x78 (27890)

Load capacity: total load 80 kg, shelf 15 kg

Material: - vertical aluminium extrusions
- plastic parts in ABS
- steel sheet epoxy powder coated

Shelf: 13 kg (27895)
15 kg (27893, 27896, 27889)
17 kg (27890)
18 kg (27897) 18.5 kg (27894)

Total weight: 13 kg (27895)
15 kg (27893, 27896, 27889)
17 kg (27890)
18 kg (27897) 18.5 kg (27894)

Medical Device Class I


GIMA code	SMART CARTS	Height (cm)
27893	2 shelves 40x36 cm + base + drawer	78
27894	2 shelves 50x42 cm + base + drawer	78
27895	2 shelves 40x36 cm + base	78
27896	3 shelves 40x36 cm + base	88
27897	4 shelves 40x36 cm + base	115
27889	2 shelves 50x42 cm + base	78
27890	2 shelves 60x42 cm + base	78

ACCESSORIES FOR SMART CARTS

27873 Basket
27874 Infusion stand - fixed height*
27888 Cardio arm for ECG

*Trolley with adjustable height infusion stand available on request in 40 days

EASY AND EASY POWER CARTS



EASY CARTS AND EASY POWER CARTS
A complete line of stylish and sturdy carts with 2, 3 or 4 shelves and with power.
Adjustable height of sturdy shelves, aluminium columns and 4 twinned castors Ø 65 mm (2 with brake).
All shelves with anti-fall edges.
Designed for medical devices.
Illustrated assembly instructions.

Easy power carts

Easy power carts are equipped with cable pipe inside lateral columns and 4 connection sockets.

TECHNICAL SPECIFICATIONS

Dimensions (WxDxH): 50x45x75/83.5/107.5 cm (Easy)
50x45x75/83.5 cm (Easy power)

Load capacity: total load 140 kg
shelf 25 kg

Shelves: 40x36 cm

Castors: 4 twinned castors Ø 65 mm (2 with brake)

Material: - vertical aluminium extrusions
- plastic parts in ABS
- steel sheet epoxy powder coated

Shelf: 13 kg (27880)
15 kg (27881, 27882, 27883)
17 kg (27886)
18 kg (27887) 18.5 kg (27884)

Medical Device Class I

GIMA code	EASY CARTS	Height (cm)
27880	2 shelves 40x36 cm + base	75
27881	3 shelves 40x36 cm + base	83.5
27882	4 shelves 40x36 cm + base	107.5

GIMA code	EASY POWER CARTS	Height (cm)
27883	2 shelves 40x36 cm + base	75
27886	3 shelves 40x36 cm + base	83.5

ACCESSORIES EASY AND EASY POWER CARTS

27873 Basket 30x16xh 14.5 cm
27874 Infusion stand - fixed height*
27884 Extra shelf 40x36 cm
27888 Cardio arm for ECG

*Trolley with adjustable height infusion stand available on request in 40 days