


SUPER EASY MULTIFUNCTIONAL CARTS

BUILD UP YOUR CART
YOU CAN ADD: EXTRA SHELF, BASKET,
INFUSION STAND, CARDIO ARM



SUPER EASY CART: A new stylish and quality cart you can build according to your needs. Carts with base and 1 or 2 sturdy shelves, aluminium columns and 4

| GIMA code | SUPER EASY CARTS - with handle | Shelf size (cm) |
|-----------|--|-----------------|
| 27870 | Cart + shelf + base | 32x26 |
| 27871 | Cart + shelf + base | 40x36 |
| 27869 | Cart + shelf + base + monitor support* | 40x36 |

*Can not be used with 27874

twinned castors Ø 65 mm (2 with brake). Illustrated assembly instructions.

| GIMA code | ACCESSORIES for SUPER EASY CARTS |
|-----------|---|
| 27872 | Extra shelf 32x26 cm |
| 27873 | Basket |
| 27874 | Infusion stand - fixed height** |
| 27875 | Extra shelf 40x36 cm |
| 27888 | Cardio arm for ECG (not suitable for 27869) |

**Trolley with adjustable height infusion stand available on request in 40 days


TECHNICAL SPECIFICATIONS

Dimensions: 47x42x79 cm (27870-1)
47x42x102 cm (27869)
Load capacity: total 120 kg - shelf 10 kg
4 twinned castors Ø 65 mm (2 with brake)
Material: - extrusions in aluminium
- plastic parts in ABS
Shelf: steel sheet epoxy powder coated
Medical Device Class I

TROLLEY FOR PATIENT MONITORS, VENTILATORS, ...

• **45760 TROLLEY FOR MONITORS - without plate adapter**
Trolley specially designed for support and mobility of patient monitors and ventilators of most manufacturers. Trolley needs plate adapter (not included) to be chosen depending on brand/model.

- 10 mm thickness steel base with high resistance ABS cover
 - column in extruded aluminium with epoxy powders oven coated surface treatment
 - ergonomic rear handle in ABS with reinforced structure
 - 5 twin pivoting anti-noise castors Ø 65 mm made of thermoplastic rubber and ABS, 3 with brake
- Storage basket made of ABS plastic (30.5x15.6xh 14.9 cm).
Size: 52.5x50xh 105 cm.
Max load: 15 kg.



| GIMA Code | Type | PLATE ADAPTERS |
|-----------|------|------------------------|
| 45761 | A | Mindray adapter |
| 45762 | B | Philips adapter |
| 45763 | C | Philips "Horn" adapter |
| 45764 | D | Other brands adapter |

| GIMA Code | Type | BRAND | MODEL |
|-----------|------|------------------|--|
| 45764 | D | Colin | Press-Mate 8800 |
| 45764 | D | Fukuda Denshi | Dynascope DS-8007 and DSA-81 |
| 45764 | D | General Electric | GE Dash 2000, 3000, 4000 monitors Carescape B650, B450 Carescape B20, B40, B105, B125 Root M series Rad-97 M series |
| 45764 | D | Masimo | ePM Compact: ePM10-12-15 ePM Modular: ePM10M, ePM12M Vet |
| 45761 | A | Mindray | iPM: iPM 8-10-12, iMEC: iMEC 8-10-12 uMEC: uMEC 8-10-12 Benevision: N12-15-17 VS600, VS900 Accutor 3-7 |
| 45764 | D | Nellcor | N-595, N-550 Pulse Oximeter Bedside SpO ₂ PM100N |
| 45764 | D | Nihon Kohden | Life Scope i, Life Scope L, BSM 2301, BSM 2303 Vismo PVM 2700, Vismo PVM 4000 Life Scope VS, BSM 3500, BSM 3700, Vismo PVM 2700 Life Scope P, BSM 4103, BSM 4113 Life Scope A, BSM 5105, BSM 5135 Life Scope TR, BSM 6301K, BSM 6501K, BSM 6701K Life Scope G5, Life Scope N, OPV 1500 GF-210R Multigas Unit CSM-1501, CSM-1502, BSM-6XXX, WEP-5200, GF-200 series Capnostream™ 20 and Capnostream™ 20p Efficia CM10-12-100-120-150 EarlyVue VS30 IntelliVue MP5-20-30-40-50-60-70 IntelliVue MX400-450-500-550-600-700-800 |
| 45764 | D | Oridion | SureSigns Vsi, VS2+ |
| 45764 | D | Philips | SureSigns VM4-6-8 and VS3-VS4-VSV series |
| 45763 | C | Philips | MIST and UltraMIST series |
| 45762 | B | Philips | S and M series |
| 45764 | D | Sanuwave | |
| 45764 | D | Sonosite | |