

DOUBLE FACE PHARMACY TROLLEY



27499

MADE IN ITALY



45730



45731

• **27499 DOUBLE FACE PHARMACY TROLLEY** - 3 large, 30 small drawers
Two trolleys in one: 3 large drawers on one side, 30 (6+9+15) monodose cassettes on other side; with locking system.

Size of cassettes:

- 6 cassettes
115x290xh 75 mm

- 9 cassettes
115x290xh 45 mm

- 15 cassettes
70x290xh 45 mm

Size of trolley:
920x590xh 810 mm

Provided with 4 revolving castors Ø 125 mm (3 with brake, 1 of trolley, antistatic).



• **45730 DOUBLE FACE PHARMACY TROLLEY** - 3 large, 22 small drawers
Two trolleys in one: 3 large drawers on one side, 22 (6+6+10) monodose cassettes on other side; with locking system.

Size of cassettes:

- 6 cassettes
115x290xh 75 mm

- 10 cassettes
115x290xh 45 mm

Size of trolley:
920x590xh 810 mm
Provided with 4 revolving castors Ø 125 mm

(3 with brake, 1 of trolley, antistatic).



• **45731 DOUBLE FACE PHARMACY TROLLEY** - 62 small drawers
Two trolleys in one: 31 (6+10+15) monodose cassettes on both sides; with locking system.

Size of cassettes:

- 6 cassettes
115x290xh 75 mm

- 10 cassettes
115x290xh 45 mm

- 15 cassettes
70x290xh 45 mm

Size of trolley:
920x590xh 810 mm

Provided with 4 revolving castors Ø 125 mm (3 with brake, 1 of trolley, antistatic).



EXCEL TROLLEY

Multilingual instructions:
GB, FR, IT, ES,
PT, DE, GR.



• **27440 EXCEL TROLLEY** - 3 shelves (60x45cm) Size: 67x69xh 80 cm.

• **27441 EXCEL TROLLEY** - 4 shelves (60x45cm) Size: 67x69xh 110 cm.
Maximum load: 80kg
Delivered in kit form.

TRIS TROLLEY



• **27438 TRIS TROLLEY**

A practical trolley for ECG, Electrosurgical unit or any other device, with 4 nylon castors and 2 lateral handles for easy moving.

Size: 49x37xh 77 cm. Shelf size: 40x30 cm.
Total load: 39 kg, max per shelf 13 kg
Colour: RAL 9002

Multilingual manual: GB, FR, IT, ES, PT, DE.
Delivered in kit form.

PROFESSIONAL CARTS



• **27885 PRO CART** 2 shelves - h 80 cm
• **27887 PRO CART** 2 shelves + keyboard and monitor support

• **27873 BASKET**

• **27874 INFUSION STAND**

TECHNICAL SPECIFICATIONS

Load capacity: - total load 240kg - shelf 35kg
Base: 60.8x50 cm
Shelves: 2 shelves 50x50 cm + base
Material: - vertical extrusions in aluminium
- plastic parts in ABS