

## DRUGS & MEDICINE TROLLEYS



27495



27497




27494

- 27494 PHARMACY TROLLEY - small - 67x64xh 100 cm with 15 3-compartment partition
- 27495 PHARMACY-TROLLEY - standard - 82x64xh 100 cm with 20 3-compartment partition
- 27497 PHARMACY TROLLEY - large - 115x64xh 100 cm with 30 3-compartment partition

A portable pharmacy, allowing to manage daily or weekly treatment of patients, with 15/20/30 3 compartments small drawers for easy drug storage (45-60-90 spaces). Painted steel frame, technopolimer surface and panels. Can be completely closed by safety key lock. Provided with 4 revolving castors Ø 125 mm (2 with brake).

- 27498 PARTITION - 3 compartments - spare

MADE IN ITALY



4 x 27498

## DRUGS & MEDICINE CABINETS



27942, 27944



27946

Compartment for medicine - spare

MADE IN ITALY



27941



**MEDICINE CABINETS 100x60xh 195 cm**  
**With drugs compartment, 60 medicine compartments and 5 adjustable shelves**  
 - made of bilaminated board (27944/5) or melamine board (27942) 20 mm thick  
 - epoxy powder painted steel base frame with 4 adjustable feet Ø 50 mm  
 - 2 external lockable doors turning through 180° with 10 medicine compartments each  
 - perimetrical anti-dust gaskets  
 - inside 2 containers with 10 medicine compartments each side turning through 90°  
 - dangerous drugs compartment 40x22x40 cm provided with lock

- 27946 COMPARTMENT FOR MEDICINE for 27941-5 - spare for cabinet sold since 2022



- 27941 COMBO MEDICINE CABINETS 96x66xh 195 cm with 30 medicine compartments one side and 6 ISO trays 1 shelf and 3 drawers other side
  - made of bilaminated board 20 mm thick
  - epoxy powder painted steel base frame with 4 adjustable feet Ø 50 mm
  - Right side**
  - 1 external lockable door turning through 180° with 10 medicine compartments each
  - perimetrical anti-dust gaskets
  - inside 1 container with 10 medicine compartments each side, turning through 90°, provided with lock
  - Left side**
  - glass door, internal side, with runners for ISO trays
  - 3 drawers on the bottom
- Supplied with 6 ISO trays different heights (2xh 5 cm, 3xh 10 cm and 1xh 20 cm) and 1 shelf.

GIMA code	MEDICINE CABINET	Material	Colour
27942	Medicine cabinet	melamine particle board	white
27944	Medicine cabinet	bi-laminated board	white
27945	Medicine cabinet	bi-laminated board	any colour*
Colours:	412	431	441
	482	484	499
	854	861	1818

\*You have to indicate requested colour with order. Available in 30 days