


## BAG HOLDER TROLLEYS



45904+45910

MADE IN ITALY



45905+45910

MADE IN ITALY



45906+45910

MADE IN ITALY

- **45904 BAG HOLDER TROLLEY** - 1 bag
  - **45905 BAG HOLDER TROLLEY** - 2 bags
  - **45906 BAG HOLDER TROLLEY** - 3 bags
- Bag holder trolleys with pedal operated opening. Rigid, sturdy, self-supporting ABS plastic base with edges. The bag holding ring (35x25.5 cm) is made of glass reinforced nylon with a rubber gasket to hold the bag in place. Two vertical tubes made of AISI 304 stainless steel are welded to the base and carry the bag holding ring.

Non marking grey rubber swivelling wheels Ø 80 mm (two with pedal operated brake).

Pedal system equipped with shock absorber for silent lid closure.

- **45910 PUSH HANDLE**
  - **27466 WRAPPING BAG** - yellow strip
  - **27462 WRAPPING BAG** - blue strip
  - **45908 WRAPPING BAG** - green strip
- White synthetic bag with strip, washable up to 200 times. Size 57xh 103 cm.



45905+27462+27466+45910

## LAUNDRY TROLLEY



45913

Instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, CZ, HU, SE, FI, HR, GR, BG, LT.

- **45913 LAUNDRY TROLLEY**
- Laundry trolley with stainless steel structure. Canvas bag and three shelves. Silent and antistatic castors Ø 7.5 cm with brake assembled on the vertical column line for better stability. Size 900x450xh 1050 mm. Supplied with canvas bag.

- **45914 CANVAS BAG** - spare



Noiseless castors with brake



Velcro for canvas bag


## SHUT BAG HOLDER



45900

MADE IN ITALY

- **45900 SHUT BAG HOLDER** - 70 l
- Trolley suitable for carrying one 70 l bag (not included). Stainless steel AISI 304 structure. Equipped with swivelling castors with brake and bumpers. For use with standard 70 l plastic bag. Can use also wrapping bag (27462, 27466, 45908)



## LINEN TROLLEYS



27467

- **27467 SOILED LINEN TROLLEY**
  - **27468 DOUBLE SOILED LINEN TROLLEY**
- Chrome plated frame with 3 or 4 castors Ø 80 mm. Removable bag. Size: Ø 46xh 90 cm. Delivered in kit form.



27468

Size:  
92x90 cm