

## LAUNDRY TROLLEYS



MADE IN ITALY

- 27464 LAUNDRY TROLLEY - painted steel
- 27465 LAUNDRY TROLLEY - stainless steel

Trolleys to collect dirty laundry while distributing clean linen to the patient. Stainless steel or painted steel frame. Equipped with 4 Ø 125 mm castors (2 with brake) and 4 rubber bumpers at corners. Size: 135x65xh 102 cm.

- 27463 SIDE PANEL - stainless steel
- 27466 WRAPPING BAG - yellow strip
- 27462 WRAPPING BAG - blue strip

washable up to 200 times, 57xh 103 cm



MADE IN ITALY

- 45925 CLOSED LAUNDRY TROLLEY - painted steel
- 45926 CLOSED LAUNDRY TROLLEY - stainless steel

Trolleys for collection of dirty laundry and distribution of clean linen to the patient. Stainless steel or painted steel frame.

Equipped with 4 Ø 125 mm anti-thread swivelling castors (2 with brake), 2 rectangular supports with rubber bag stop ring, 2 lids and 3 shelves of which 2 adjustable.

Suitable for 70 litres waste bags.  
Size: 1330x650xh 1200 mm (45925)  
1330x650xh 1140 mm (45926)


- 27466 WRAPPING BAG - yellow strip
  - 27462 WRAPPING BAG - blue strip
- White washable laundry synthetic bag with yellow or blue strip.  
Size: 57xh 103 cm.

MADE IN ITALY

- 45915 LAUNDRY TROLLEY - laminated
- Closed trolley with doors, 3 shelves plus top and 2 bag holders. Laminated particle board body and stainless steel AISI 304 structure. Equipped with swivelling castors, 2 with brake and bumpers.

Size: 110x66x132 cm.

- 27462 WRAPPING BAG - blue strip
  - 27466 WRAPPING BAG - yellow strip
- White washable laundry synthetic bag with yellow or blue strip.  
Size: 57xh 103 cm.



## CLEAN LINEN AND PATIENT HYGIENE TROLLEYS



### CLEAN LINEN TROLLEYS 2 DOORS

Solid and sturdy trolleys for distribution of clean linen and patient hygiene. Steel painted with antimicrobial epoxy powders to obtain maximum lightness. Upper part has large support surface in technopolymer plastic material. Integrated ergonomic push handle. Supplied with 2 doors with 3 internal shelves (45930, 45933 only) or 2 doors with 2 internal shelves and one drawer made of antimicrobial epoxy powder coated steel (45931, 45934 only). Key closures.

4 castors Ø 125 mm, anti-wire, pivoting, two of which with lever locking system and reinforced plastic support.

GIMA code	CLEAN LINEN TROLLEYS	Size (mm)	Drawer size (mm)
45930	Clean linen trolley - small	670x640xh 1280	-
45931	Clean linen trolley - small	670x640xh 1280	600x480xh 150
45933	Clean linen trolley - standard	820x640xh 1280	-
45934	Clean linen trolley - standard	820x640xh 1280	750x480xh 150