


ELECTRIC CHAIRS WITH FOOT CONTROL - 4 MOTORS



Size: 185x88xh 63-90 cm
Net weight: 98 kg
Max load: 175 kg

- 28047 SABA ELECTRIC CHAIR - 4 motors - white
 - 28048 SABA ELECTRIC CHAIR - 4 motors - avio blue
- Electrically controlled 4 motors chair to adjust separately angle of the backrest, seat inclination, footrest and lifting of the whole chair. Hand controller and built-in foot controller which allow operator to adjust chair. High-grade upholstery leather easy to clean and fresh feeling. Headrest with face hole. Supplied with removable cushion and remote controller. Metal base for high stability. Operating power: 100-240V - 50/60Hz. Manual: GB, FR, IT, ES.

- 28013 COUCH ROLL HOLDER - optional
- 45231 STOOL - white
- 45235 STOOL - avio blue
- 45241 STOOL WITH BACKREST - white
- 45245 STOOL WITH BACKREST - avio blue



28048

PODOLOGY CHAIRS

Size: 186x86xh 58-89 cm
Net/Gross weight: 72/92 kg
Max load: 175 kg
Operating power: 220-240 V - 50/60 Hz.



45232

28025

- 28025 ELECTRIC PODOLOGY CHAIR - 3 motors - blue
- Electrically controlled chair to independently adjust height, backrest and seat inclination, by optional foot switch or push button panel. Headrest can be inclined, arm rests are movable. Manual: GB, IT.
- 28013 COUCH ROLL HOLDER
 - 45232 STOOL - blue
 - 45242 STOOL WITH BACKREST - blue

Size: 180x80xh 51-65 cm
Net weight: 41 kg
Max load: 175 kg



45241


28020

- 28020 PODOLOGY CHAIR - mechanical - white
- This chair can move up and down and rotate. Backrest and footrests are adjustable. Heavy round base provides maximum stability. Manual: GB, IT.
- 45231 STOOL - white
 - 45241 STOOL WITH BACKREST - white

PODOLOGICAL CABINETS

- 28030 GAMMA1 CABINET - 61x45xh 77 cm - small
 - 28031 GAMMA3 CABINET - 108x45xh 77 cm - large
- Equipped with tabletop in scratch-resistant wood, drawers and doors, designed to keep tools sterile. Supplied with connections to the mains: germicide compartment, lampholder and other electric sockets. Power supply: 220 V - 50/60 Hz.

	Cabinet		Table top		Germicide compartment	Large compartment
	28030	28031	28030	28031		
Length:	61 cm	108 cm	61 cm	108 cm	28 cm	28 cm
Width:	45 cm	45 cm	42 cm	42 cm	38 cm	28 cm
Height:	77 cm	77 cm			4 cm	15 cm



28030

28031