

GIMA ADJUSTABLE HEIGHT ELECTRIC TABLES - FOOTSWITCH

Sturdy and highly stable examination and treatment table with white metal painted frame. Adjustable height by electric footswitch.

Fixed leg-section and manually adjustable head-section with mouth-nose hole.

Available in 12 colours. Fireproof padding.

Class I Medical Devices, complies with MDR 2017/745.

Lifting capacity: 135 kg. Length: 193 cm. Weight: 54 kg.

Power supply 100-240 V, 50-60 Hz.


Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Matching stools see page 181.

ADJUST. HEIGHT ELECTRIC TABLE - FOOTSWITCH				
Electric table		Colour	No.	Mouth-nose hole plug - optional
27628	Electric table - footswitch	green	4	27625
27633	Electric table - footswitch	blue	2	27631
27634	Electric table - footswitch	beige	5	27620
27635	Electric table - footswitch	black	1	27621
27632	Electric table - footswitch	any colour*	other	27630

*Available in 20 days, see beside available colours

ACCESSORIES	
27629	Couch roll holder
27650	Foot Switch - spare



GIMA ADJUSTABLE HEIGHT ELECTRIC TABLES - FOOTBAR

Sturdy and highly stable examination and treatment table.

Adjustable height by electric footbar.

Adjustable head-section by gas pump with mouth-nose hole.

Sturdy tubular steel frame painted with non-toxic epoxy powder, RAL 9010. Chromed footbar. Intended for professional use in hospitals, physiotherapy, treatments. Supplied with 6 cm thickness mattress and mouth-nose hole plug 23x11 cm included.

Fireproof padding. Available in 12 colours.

Class I, Medical Devices, complies with MDR 2017/745.

Maximum overall dimensions 195x59,5xh 54-71 cm.

Weight: 65 kg. Lifting capacity: 135 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

ADJUST. HEIGHT ELECTRIC TABLE - FOOTBAR			
Electric table		Colour	No.
44520	Electric table - footbar	blue	2
44521	Electric table - footbar	beige	5
44522	Electric table - footbar	black	1
44525	Electric table - footbar	any colour*	other

*Available in 20 days, see beside available colours

ACCESSORIES	
27629	Couch roll holder



VARIABLE HEIGHT EXAMINATION COUCH WITH TR/RTR

• 44510 LORD VARIABLE HEIGHT EXAMINATION COUCH - water green - 5

• 44516 LORD VARIABLE HEIGHT EXAMINATION COUCH - any colour*

*available in 20 days, see below available colours

Stable and sturdy examination couch made of high quality tubular for furniture. Very practical and user friendly. Height adjustable thanks to the two electrical columns that can move simultaneously or separately (allowing the Trendelenburg/Reverse Trendelenburg position)

by a push-button panel. Under the top frame, an ergonomic lever activates a gas spring that reclines the cushions. High density, rubber padded, fire retardant and Class 1 IM top, covered in coated fabric (vinyl, synthetic leather, flame retardant Class 1 IM). It can be cleaned with common detergents.

Complies with: EU 93/42 CEE.

Maximum load: 200 kg.

Supplied with mouth/nose hole plug, paper roll holder and 4 adjustable feet.

• 44505 RETRACTABLE/REMOVABLE/TURNABLE WHEELS KIT - optional

