


TEST TUBES, MICROTUBES AND PIPETTES

In order to fulfill every laboratory operator's necessity, a wide range of disposable test tubes, microtubes and pipettes is available in different sizes and materials. High quality raw materials, transparency and resistance are the main technical features that characterize these products.



- maximum transparency
- excellent resistance
- accuracy of graduation
- wide range of products

POLYSTYRENE TEST TUBES WITH OR WITHOUT RIM



Stoppers in polyethylene, neutral colour



GIMA code	POLYSTYRENE TEST TUBES	Minimum order
TEST TUBES WITHOUT RIM		
37902	Test tube 5 ml - cylindrical - 12x75 mm	box of 4,000*
37903	Test tube 5 ml - cylindrical - 13x75 mm	box of 4,000*
*16 bags of 250 pcs		
GRADUATED TEST TUBES WITH RIM		
37912	Test tube 10 ml - cylindrical - 16x100 mm	box of 2,000
37913	Test tube 10 ml - conical - 16x100 mm	box of 2,000
WINGED STOPPERS FOR TEST TUBES		
37930	Stopper for tubes Ø 12 mm	box of 1,000
37931	Stopper for tubes Ø 13 mm	box of 1,000
37932	Stopper for tubes Ø 16 mm	box of 1,000

STERILE TEST TUBES - WITH RED CAP



Versatile test tubes with red caps, which indicate that tubes are sterile.


All tubes are provided with writing area.

Beta rays sterilization.

GIMA code	STERILE TEST TUBES	Minimum order
POLYPROPYLENE TEST TUBES		
37920	Test tube 5 ml - cylindrical - 12x86 mm	box of 1,500
37921	Test tube 10 ml - cylindrical - 16x100 mm	box of 800
POLYSTYRENE TEST TUBES		
37925	Test tube 5 ml - cylindrical - 12x86 mm	box of 1,500
37926	Test tube 10 ml - cylindrical - 16x100 mm	box of 800

MICROTUBES

Graduated microtubes made of polypropylene. Flat cap (37941-2 only), one handed usage, neutral colour, frosted area and perfect closure.



GIMA code	CONICAL MICROTUBES	Minimum order
37940	Microtube 0.5 ml	box of 16,000*
37941	Microtube 1.5 ml	box of 10,000**
37942	Microtube 2 ml	box of 8,000***


*16 pouches of 1000 ** 20 pouches of 500

***16 pouches of 500

TEST TUBE RACKS

- 37935 RACK FOR TEST TUBES Ø 12 AND 13 mm - 90 places - box of 5

Autoclavable, made of polypropylene



- 37936 RACK FOR TEST TUBES Ø 16 mm - 60 places - box of 5


Autoclavable, made of polypropylene



PASTEUR PIPETTES



Pipette 1 ml 37950, 37953



Pipette 3 ml 37951, 37954



37953, 37954


Made of polyethylene, individually wrapped in peel pack

GIMA code	PASTEUR PIPETTE	Minimum order
37950	Pipette 1 ml	box of 5,000*
37951	Pipette 3 ml	box of 5,000*
37953	Pipette 1 ml - sterile	box of 1,000
37954	Pipette 3 ml - sterile	box of 1,000

*10 boxes of 500 pieces

CENTRIFUGES TEST TUBES - GRADUATED WITH CAP

All test tubes have a writing area and screw cap



GIMA code	CENTRIFUGES TEST TUBES IN POLYPROPYLENE	Minimum order
TEST TUBES - NON STERILE		
37960	Test tube 15 ml - conical - 17x120 mm	box of 1,700
37961	Test tube 50 ml - conical - 30x115 mm	box of 500
37962	Test tube 50 ml - conical with base*	box of 500
TEST TUBES - STERILE		
37965	Test tube 15 ml - conical - 17x120 mm	box of 1,700
37966	Test tube 50 ml - conical - 30x115 mm	box of 500
37967	Test tube 50 ml - conical with base*	box of 500

*Size 30x115 mm