

WING SCREEN WITH TREVIRA CURTAINS



SCREEN WITH TREVIRA CURTAINS

Build up the screen you need choosing among 4 frames and 6 colours of Trevira curtains.

Wing screen

Aluminium alloy tubular frame mounted on feet or castors (Ø 50 mm).

Size: 150xh 167.5 cm (45560)

150xh 174 cm (45561)

200xh 167.5 cm (4 wings)


Packing size: 176x57x13 cm (all wing screens)

Trevira curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C.

Available in 6 colours to match any need.

Custom-made curtains available on request.

Multilingual manual: GB, FR, IT.



GIMA Code	WING SCREEN FRAME without curtain
45560	3 wing screen
45561	3 wing screen with castors
45562	4 wing screen
45563	4 wing screen with castors
GIMA code	CURTAINS* IN TREVIRA 45xh 129 cm
45570	Curtains - light blue
45571	Curtains - blue
45572	Curtains - white
45573	Curtains - beige
45574	Curtains - green
45575	Curtains - peach


*you have to buy 3 or 4 curtains depending on frame you choose



MADE IN ITALY



3 OR 4 WINGS SCREEN WITH PLASTIC CURTAINS



27970



27974

- **27970 3 WINGS SCREEN**
 - **27971 3 WINGS SCREEN with castors**
 - **27972 4 WINGS SCREEN**
 - **27974 4 WINGS SCREEN with castors**
- Aluminium alloy tubular frame.
Weight 5 kg (3 wings), 8 kg (4 wings).
Foldable washable white plastic curtains mounted on feet or castors (Ø 50 mm).
Size: 150xh 170 cm - 3 wings (50 cm each)
200xh 170 cm - 4 wings (50 cm each)
- **27975 KIT OF 6 CASTORS - spare for 27971/4, optional for 27970/2**

MADE IN ITALY

SEPARATION SYSTEMS AND CURTAINS see page 139

2 OR 3 WINGS WARD SCREENS WITH TREVIRA OR COTTON CURTAINS



45580

Height: 180 cm



45579

Width: 60+60 cm

- **45578 2 WINGS SCREEN - white cotton cloth**
 - **45579 2 WINGS SCREEN - light blue Trevira**
 - **45580 3 WINGS SCREEN - white cotton cloth**
 - **45581 3 WINGS SCREEN - light blue Trevira**
- Exclusive 2 or 3 wings ward screens with sturdy structure in oval tube (30x15 mm), painted with epoxidic powder metal grey.
Light blue Trevira curtains or white washable cotton cloth.
Equipped with 4 adjustable feet that ensure increased stability.
Class I Medical Device complying with European Directive 93/42

MADE IN ITALY

Width: 60+60+60 cm