

## CABINETS FOR MEDICINES AND INSTRUMENTS

CABINETS ARE  
SUPPLIED  
DISASSEMBLED

MADE IN ITALY



**27898**

Size: 52x35x138 cm  
Pack size: 65x26x134 cm  
Weight: 31 kg



**27899**

Size: 75x37x164 cm  
Pack size: 72x24x156 cm  
Weight: 45 kg



**27900**

Size: 75x37x164 cm  
Pack size: 72x24x156 cm  
Weight: 50 kg



**27902**

Size: 79x36x186 cm  
Pack size: 88x36x96 cm  
Weight: 56 kg

- **27898 VALUE CABINET - 1 door**
- **27899 VALUE CABINET - 2 doors**

Cabinets made of enamelled steel sheet with tempered glass doors. Hinged glass door with lock and key. Three stainless steel shelves adjustable in height. **Delivered in kit form.**

- **27900 VALUE CABINET - 4 doors (hinged)**
- **27902 VALUE CABINET - 4 doors (2 hinged, 2 sliding)**

Cabinets made of enamelled steel sheet with 2 overlapped compartments. Upper compartment with 2 tempered glass doors (hinged doors 27900 or sliding doors 27902), lock with key and 2 stainless steel, height adjustable shelves. Inferior compartment with 2 doors made of enamelled steel plate and 1 enamelled steel sheet shelf. **Delivered in kit form.**

## CABINETS WITH 1, 2, 4 DOORS AND WALL CABINETS

MADE IN ITALY



**27901**

**27905**

**27903**

Lower section  
with 1 steel shelf



**27912**

Hanging cabinet with 2 lockable doors. Inside division in 2 compartments left side with 2 removable shelves, right side with 1 removable shelf.



**27911**

Hanging cabinet to keep medicine with 2 glass tempered doors and 2 tempered glass shelves adjustable on master rack to be hung by two holes on the back. Supplied with lock.

**GIMA CABINETS** Choice of cabinets with epoxy coated steel structure and tempered glass. Rack-mounted adjustable glass shelves. Rubber feet.

TECHNICAL SPECIFICATIONS	27901	27905	27903	27911	27912
Door:	1	2	4	2	2
Adjustable glass shelves:	3	4	2	2	3
Weight kg:	30	47	45	16	18
Size cm:	53x36x144	75x38x169	75x38x169	60x30x60	80x33x45

- **27901 CABINET - 1 door**
- **27905 CABINET - 2 doors**
- **27903 CABINET - 4 doors**
- **27911 WALL CABINET - glass sliding door**
- **27912 WALL CABINET - painted door**