

OP-DOP VISOR SHIELD "PROTECTOR" WITH 3+1 SHIELDS

• 25642 VISOR SHIELD "PROTECTOR" - white

The ultra-light face shield that combines function with style:

protection: a waterproof plastic barrier will protect you

comfort: fast and easy to wear and remove

convenience: reusable and replaceable shield

optical quality: superior optical quality.

Medium adjustable size: fits over 90% of healthcare professionals. Fits over glasses.

Shielded air vent eliminates fogging.


Supplied with 3 face transparent shields + 1 head shield.

• 25643 REPLACEMENT SHIELD KIT

Includes 2 face transparent shields

+ 1 head shield + 1 foam bumper.

MADE IN THE USA



25643



CE

25642

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, GR.



VISOR HEAD MASKS

VISOR HEAD MASKS

Great features

- faceshield with 90° rotation
- inner and outer surface treated against mist
- superior coverage to protect against droplets and splashes

Superior materials

- elastic band in anallergic material
- rotation systems in plastic material
- soft support in anallergic interchangeable foam
- faceshield in heat-formed acetate with anti-fog treatment

Easy and comfortable use

- fitting headbands adjustable by velcro fasteners
- usable over ophthalmic glasses

GIMA code	VISOR HEAD	Minimum order
25656	VISOR HEAD MASK	1
25668	VISOR SHIELD	1
25669	SPARE SHIELD for 25668	box of 5

Antifog mask
0.5 mm thick

CE

Class II EN 166-167-168 U1SN



25656



25668

Faceshield with 90°
max rotation



Antifog mask
0.5 mm
thick

25668

GIMA PROTECTIVE SHIELDS


• 30846 GIMA PROTECTIVE SHIELD - fixed shield

• 30843 GIMA PROTECTIVE SHIELD - liftable shield

Protective visor made of polycarbonate with large transparent surface and high coverage degree, it allows to maintain eye contact by protecting eyes and mucous membranes in complete safety.

Sanitizable, universal size and extremely ergonomic thanks to its simil velcro closure band and soft, interchangeable pad.

MADE IN ITALY



30843 Weight 125 g
30846

30843
30846

CE 0426

30843
Liftable shield

Suitable for clinical/professional and for personal use.

It is usable with glasses, wide peripheral vision.

Gima protective shield is compliant with International Standards (PPE, UE 2016/425, Class II, EN 166:2001).

Supplied in a polyethylene bag with multilingual instructions: GB, FR, IT, ES, PT, DE, GR, RO, PL, SE, HU, CZ.

• 30847 FOREHEAD PADDING for 30843, 30846 - spare



30846

Interchangeable
forehead padding

DISPOSABLE ANTI-FOG FACE SHIELD

• 25662 DISPOSABLE ANTI-FOG FACE SHIELD

Protective disposable face shield made of anti-fog hydrophilic film. Clear and smooth visor allows to maintain eye contact by protecting whole face against potential contamination from blood pathogens, fluids or harmful chemical splash.

It is equipped with sponge strip and elastic band.

Size: 330x220 mm.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, RO, PL.

17



Sponge strip:
270x35x35 mm
Elastic band:
12 mm

CE

25662

