


HEIGHT ADJUSTABLE STOOLS



- 45210 STOOL - black
 - 45211 STOOL - blue
 - 45212 STOOL - white
- Padded simulated leather seat Ø 36 cm. Height adjustable 46-60 cm. Plastic base with 5 castors Ø 50 mm.
- 45220 STOOL - black
 - 45221 STOOL - white
 - 45222 STOOL - red
- Padded fire proof simulated leather seat and back. Height adjustable. Plastic base with 5 castors and chromed ring.



- 27949 STOOL - feet
 - 27954 STOOL - castors
- Stainless steel seat AISI 304 Ø 32 cm. Height adjustable 45-64 cm (27949), 50-68 cm (27954).



- 27947 STOOL - feet
- Stainless steel seat AISI 304 Ø 32 cm. Height adjustable 46-65 cm. Aluminium base with feet.



- 27955 STOOL - grey
- Washable padded simulated leather seat Ø 31 cm. Height adjustable 48-67 cm. Plastic base with feet.



- 27956 STOOL - grey
- Washable padded simulated leather seat Ø 31 cm. Height adjustable 53-71 cm. Plastic base with 5 castors Ø 50 mm.

HEIGHT ADJUSTABLE STOOLS



- 45230 STOOL - black
 - 45231 STOOL - white
 - 45232 STOOL - blue
 - 45235 STOOL - avio blue
 - 45236 STOOL - light beige
- Padded simulated leather seat Ø 34 cm. Height adjust. 55 to 75 cm.
- 45227 CASTOR PLUG-IN TYPE - spare
 - 45228 CASTOR SCREW TYPE - spare
- 45240 STOOL - black
 - 45241 STOOL - white
 - 45242 STOOL - blue
 - 45245 STOOL - avio blue
 - 45246 STOOL - light beige
- Padded simulated leather seat Ø 34 cm with backrest. Height adjustable 55 to 75 cm.
- 45227 CASTOR PLUG-IN TYPE - spare
 - 45228 CASTOR SCREW TYPE - spare

ERGONOMIC HEIGHT ADJUST. STOOLS



- 45250 ERGO STOOL - black
 - 45251 ERGO STOOL - white
 - 45252 ERGO STOOL - blue
- Ergonomically designed comfortable stool with saddle-shaped seat and height adjustment with a gas piston. Height adjustable 52 to 67 cm.
- 45227 CASTOR PLUG-IN TYPE - spare
 - 45228 CASTOR SCREW TYPE - spare
- 45260 ERGO B STOOL - black
 - 45261 ERGO B STOOL - white
 - 45262 ERGO B STOOL - blue
- Ergonomically designed comfortable stool with saddle-shaped seat. Three levers with a gas piston allow adjustment of height, backrest angle and seat. Height adjustable 63 to 83 cm.
- 45227 CASTOR PLUG-IN TYPE - spare
 - 45228 CASTOR SCREW TYPE - spare