


GOGGLES KIT



- **25646 GOGGLES KIT - fog resistant**
Practical kit composed of 5 supports in different colours and 20 antifog shields.
Usable over ophthalmic glasses.
Norms: EN 166, EN 167, EN 168.
- **25647 SPARE SHIELDS - box of 20**

AVATAR GOGGLES



- **25264 AVATAR GOGGLES - black/red - fog resistant, antiscratch**
Anti-fog coating goggles provide anti-scratch performances to keep lenses clear and workers safe. Ultra-light and soft brow frame. 9 base lenses made of polycarbonate with hydroshield coating. Provided with innovative and adjustable nose pad.
PPE Category II. Norms: EN 166, EN 170.

HIGH PROTECTION GOGGLES AND MASK



- **25266 MEDICAL ISOLATION GOGGLES - fog resistant and anti-scratch - disposable - box of 10**
Fully enclosed goggles with air vents for protection against: droplets, splashes of liquids and body fluids. Wrap around eyewear with side shields for added protection. Clear constructed lenses for unobstructed sight. Made of transparent PVC frame and PC ocular, white PP air vents. Provided with elastic band. Non-sterile. Conform with PPE regulation (EU) 2016/425 and EN 166:2001. Cat. II. Individually packed with multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, NO, PL, CZ, SI, HR, HU, RO, BG, GR, Arabic.



- **25265 HIGH PROTECTION MASK 6x1 - fog resistant, anti-scratch**
 - wide lens for a field of view without distortion
 - temples adjustable in length and inclination with SoftPad technology
 - sophisticated indirect ventilation system
 - soft internal rubber for comfort without local pressure
 - certified according to EN166 point 3 (droplets and splashes) and point 4 (large dust particles)
 - supplied with elastic band

MASK-OP® 99% FILTERING MASKS WITH VISOR - TYPE IIR

MASK-OP® WITH VISOR

- **25652 MASK-OP® - ear loops - green - box of 50**
- **25654 MASK-OP® - lacets - light blue - box of 50**

MASK-OP® WITH WRAPPING VISOR

- **25636 MASK-OP® - ear loops - green - box of 50**
- **25637 MASK-OP® - lacets - light blue - box of 50**

Liquid splash resistant and anti-fogging surgical mask with antistatic visor.

- 4 layers: three made of non-allergenic woven non-woven and one of polyethylene net
- nuisance-free and adaptable nose-clip
- filtering efficiency (B.F.E.) $\geq 99\%$, according to EN 14683:2019+AC:2019
- resistance to breath $<30,0 \text{ Pa/cm}^2$, according to EN 14683:2019+AC:2019
- latex-free
- non sterile
- liquid splash resistance: resistant at pressure of 16 kPa, according to ISO 22609, related to EN 14683:2019+AC:2019
- anti-fogging internal stripe made of expanded polyethylene
- Bioburden $\leq 30 \text{ CFU/g}$, according to EN 14683:2019+AC:2019

Medical Device Class I Regulation EU 2017/745.

Multilingual box: GB, FR, IT, ES, PT, DE, GR, PL, RU, Arabic.

MASK-OP®
WITH
VISOR



25652

B.F.E. $\geq 99\%$

CE

MADE IN ITALY

MASK-OP® WITH
WRAPPING
VISOR



25637



25654



25636