

## GIMA FOETAL MONITOR PAPERS FOR DIFFERENT BRANDS

GIMA Code	FOETAL MONITOR manufacturer and model	Standard FHR rate/BPM**	Type	Width mm	Length mm or m'	Sheet	Grid***	Min. order
<b>BIONET</b>								
29526	Fetalcare FC/1400	Int 60/200	roll	152	25	-	green	10
29522	Fetalcare FC/700	US 60/240	roll	215	23	-	green	5
<b>COMEN</b>								
29556	Star 5000/C	Int 60/210	pack	112	100	150	or/red	20
29557	Star 5000/D	US 30/240	pack	151	100	150	or/red	25
<b>COROMETRICS</b>								
29561	4305-BAO mod. 116-145-150-151-118	Int 60/200	pack	152	90	150	or/red	25
29562	4305-BAO mod. 116-145-150-151-118	US 30/240	pack	152	90	150	or/red	25

\*mm (pack) or m (roll)

\*\*International or US standard FHR rate/BPM


GIMA Code	FOETAL MONITOR manufacturer and model	Standard FHR rate/BPM**	Type	Width mm	Length mm or m'	Sheet	Grid***	Min. order
<b>EDAN</b>								
29569	Cadence II	US 30/240	pack	112	90	150	or/red	20
29570	F/6-F/9	Int 60/200	pack	152	90	150	or/red	25
29571	Cadence II	Int 60/200	pack	112	90	150	or/red	20
<b>HEWLETT PACKARD</b>								
29576	M1911A	Int 60/200	pack	151	100	150	green	25
<b>SHENZHEN CREATIVE</b>								
29556	PC-8000 PRO single	Int 60/210	pack	112	100	150	or/red	20
29558	PC-8000 PRO twins	Int 60/210	pack	151	100	150	or/red	25
29557	PC-8000 PRO twins	US 30/240	pack	151	100	150	or/red	25
<b>SONICAID</b>								
29580	Meridian/800 - Team	White paper	pack	143	150	300	white	10

\*\*\*or/red: grid could be orange or red colour

## UNIVERSAL TRANSDUCER BELTS FOR FOETAL MONITORS

### • 29544 SINGLE PATIENT FOETAL BELTS - 100 pairs (1 pink +1 light blue)

Single patient elastic foetal belts to avoid cross-contamination. Soft and comfortable belts resistant to abrasion, corrosion, light and heat. Great elasticity and insulation. Performance: more than 1,000 stretches. Easy to wash, made of 85% polyester and spandex. Supplied with buttons.



### • 29545 REUSABLE FOETAL TRANSDUCER BELTS - box of 10

Reusable elastic fetal belt resistant to abrasion, corrosion, light and heat. Soft, comfortable and easy to wash and dry. Performance: more than 100,000 stretches.

Made of 85% polyester and spandex. Supplied with buttons.

**COMPATIBILITY (29544,29545):** Bionet, EDAN, Philips M1350, Philips FM 20 series, GE, Corometrics, Sonicaid, Huntleigh and most foetal monitors.



## OBSTETRIC STETHOSCOPES

### • 32501 SHORT OBSTETRIC STETHOSCOPE - beech-wood

### • 32502 LONG OBSTETRIC STETHOSCOPE - beech-wood

These stethoscopes are hand-craft manufactured, made of beech-wood accurately selected and seasoned for at least 3 years. They have a perfect finish with a special transparent varnish. Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, RO, HU, GR, Arabic.

### • 32505 PINARD STETHOSCOPE - aluminium



## LATEX PROBE COVERS



### • 29509 LATEX PROBE COVERS - transparent - bag of 1000 pcs - bulk

### • 29513 LATEX PROBE COVERS - transparent - box of 50

### • 29514 LATEX PROBE COVERS - transparent - box of 500

Probe covers for ecography and doppler. Length 17 cm, Ø 3.3 cm. Individually packed in single plastic bag, except 29509.

## ENEMA IRRIGATOR

### • 28720 ENEMA IRRIGATOR- 2 l capacity

Supplied in cardboard box with 2 cannulas (vaginal, anal), tap with opening lever, tube and carrying bag.

Multilingual box: GB, FR, IT, ES.

