

S/S GINAECOLOGY INSTRUMENTS AND KITS



26817

• **26817 HEGAR DILATORS SET - chrome plated**
Include 14 dilators Ø 4-17 mm




26810

• **26810 GYNAECOLOGY KIT**
Includes 4 instruments in a nylon bag. See page 45.

- 26801 FOERSTER FORCEPS 25 cm
- 26802 SCHROEDER FORCEPS 25 cm
- 26806 SCHUMACHER BIOPSY FORCEPS 24 cm
- 26808 SCHUBERT BIOPSY FORCEPS - 21 cm
- 26809 KOGAN ENDOSPECULUM - 24 cm

OTHER GYNAECOLOGY INSTRUMENTS
See page 45


S/S DELIVERY SURGICAL SET



26966

• **26966 MINOR DELIVERY SET - 25 pcs**
Kidney Dish 20 cm (8") - medium
Kidney Dish 25 cm (10") - large
Gallipot 150 ml - 6 oz - medium
Brown Sadler episiotomy scissors - 14 cm
Kocher artery forceps straight - 16 cm
Kocher artery forceps straight - 18 cm
Kocher artery forceps curved - 16 cm
English TOE dissecting forceps plain - 14 cm
English TOE dissecting forceps plain - 18 cm
Dissecting forceps toothed 1x2 - 14 cm
2 operating scissors sharp/blunt - straight - 15 cm
2 operating scissors sharp/blunt - curved - 15 cm
2 rampley sponge holding forceps - str. - 24 cm
Museum Vulsellum forceps 8 mm - str. 24 cm
3 Cusco vaginal specula - small, medium, large
3 Sims vaginal specula - small, medium, large
Umbilical scissors USA Patent - 10.5 cm
Instrument tray with lid 45x30x6.2 cm

UMBILICAL CLAMP




MADE IN EUROPE

29942

Extra safe clamp avoid accidental re-opening. Can be re-opened only by cutting back side

• **29942 X-SAFE UMBILICAL CORD CLAMP - sterile - box of 500**
Disposable umbilical cord clamp to allow a proper clamping of all umbilical cords and prevent their outcoming from the closing comb. Individually packed in coupled cellophane (nylon/polyethylene) and medical paper to ensure safe hygienic conditions to user. Sterilized by ethylene oxide.

GYNAECOLOGY SET - STERILE



Non woven bi-layer fabric gynaecology set

• **23641 GYNAECOLOGY SET - sterile - box of 8**
1 table cover 150x200 cm with central absorbent band
1 Mayo stand cover 80x145 cm
1 gynaecology drape for lithotomy position 290x240 cm with adhesive hole, built-in leggings, bag for liquid collection with filter and drain valve
1 drape 75x90 cm
1 incision film 35x50 cm
3 adhesive strips 10x50 cm
4 towels 35x40 cm

• **23644 SET FOR CESAREAN SECTION - sterile - box of 8**
1 table cover 150x200 cm with central absorbent band
1 Mayo stand cover 80x145 cm
1 drape for caesarean section
260x200x310 cm with incision film and bag for liquid collection
1 drape "baby blanket" 75x90 cm
See picture at page 40


MEMBRANE PERFORATOR



29945

• **29945 AMNIOTIC MEMBRANE PERFORATOR - sterile - box of 600**
Disposable amniotic perforator. This system is thought to be completely safe to avoid lesions to the foetus. Supplied in single sterile envelope.

I.D. BRACELET




29937

29938

GIMA code	I.D. BRACELET	Colour	Minimum order
29936	Newborn	light blue	box of 100 pcs.
29937	Newborn	pink	box of 100 pcs.
29938	Mother	light blue	box of 100 pcs.
29939	Mother	pink	box of 100 pcs.

CUPS FOR HYSTEROGRAPHY



CE 0459

29923 Connector for suction tube

29921 Cup for hysteroigraphy

• **29920 CUPS FOR HYSTEROGRAPHY - 25 mm - sterile**
• **29921 CUPS FOR HYSTEROGRAPHY - 30 mm - sterile**
• **29922 CUPS FOR HYSTEROGRAPHY - 35 mm - sterile**
• **29923 CUPS CONNECTOR - sterile**

Rigid, transparent tapered cup with a central conical cannula. Equipped with a 25 cm long flexible, pliable manipulator joined to the cup and two tubes equipped with clamps. One of these provides the vacuum, while the other injects the radiopaque medium. Perfect adhesion to the cervix by self-centering is ensured by the tapered shape of the cup and the progressive application of a vacuum to the cup. The suction cup/cannula unit is radiotransparent and allows the proximal area of the endocervix to be controlled. The radiologist is able to manipulate at a distance from the radioactive field while at the same time comfortably immobilizing the uterus. Latex free.