



# IUD - INTRAUTERINE CONTRACEPTIVE DEVICE



**• 29710 INTRAUTERINE DEVICE TCU380 PLUS - sterile**  
Sterile disposable medical device Class III for intrauterine contraception in women of child bearing age. Can be left inserted for a maximum of 5 years. Packed in a pouch with an insertion tube equipped with a movable flange to aid in gauging the right depth of insertion and a rod. Complies with ISO-7349. Multilingual manual: GB, FR, IT, ES, PT.




**• 29712 GYNEAS IUD KIT - sterile**  
A complete disposable kit for IUD insertion including: 1 intrauterine device (IUD), 1 speculum medium size, 1 hystrometer, 1 pair of Pozzi tenaculum forceps, 1 Evolution forceps, 1 pair of thread-cutting scissors, 1 practitioner and 1 patient brochure.




**• 29714 IUD INSERTION AND REMOVAL KIT - sterile**  
Disposable kit for IUD insertion and removal including: 1 Evolution Pozzi forceps, 1 long Littauer scissors, 1 Evolution long forceps.

## GYNAECOLOGY FORCEPS AND SCISSORS

**• 29717 CHERON FORCEPS 25 cm - box of 25**  
Suitable for:  
- cervix cleaning/disinfection before any intervention (ex. IUD insertion) used with a gauze or a piece of cotton and disinfectant solution  
- IUD removal




**• 29719 EVOLUTION LONG FORCEPS 25 cm - box of 25**  
Same applications as Cheron forceps, but with longer and thinner jaws.




Sterile disposable medical device, Class I. Latex free. Sterilized by ETO. Shelf lifes 5 years.

**• 29718 POZZI FORCEPS 25 cm - box of 25**  
Ergonomic, strong and thin single-use forceps suitable for grasping of the cervix, cervical cerclage, biopsies, conisations, IUD insertion (align cervix and womb with vagina axis).




**• 29720 LONG LITTAUER SCISSORS 22 cm - box of 25**  
Suitable for cutting the IUD threads after the insertion.



## BIOPSY FORCEPS

**• 29724 GYN&PUSH CERVICAL BIOPSY FORCEPS 21.5 cm - box of 25**  
Effective cervical biopsy forceps for a successfully biopsy, quickly achieved with a high quality sample. Complete comfort and safety for both doctor and patient. Disposable sterile medical device, Class II a. Latex free and phtalates free. Shelf life: 5 years.



## FLEXIBLE HYSTEROMETERS

**SINGLE USE, STERILE, FLEXIBLE HYSTEROMETERS**  
Hysterometers with soft distal part on a rigid handle for IUD\* (29925, 29926). Hysterometer with progressive dilation distal part from CH 10 to CH 14 on 25 mm, suitable for IUS\*\* (29927).



| GIMA code | HYSTEROMETERS SINGLE USE - STERILE - LATEX FREE | Ø       | Minimum order |
|-----------|---|---------|---------------|
| 29925     | Hysterometer CH 10 - 30 cm                      | 3.3 mm  | box of 25     |
| 29926     | Hysterometer CH 12 - 30 cm                      | 4.3 mm  | box of 25     |
| 29927     | Hysterometer CH 14 - 30 cm                      | 4.67 mm | box of 25     |

\*IUD: Intra Uterine Device \*\*IUS: Intra Uterine System

## EVOLUTION HYSTEROMETER

**• 29929 EVOLUTION DOUBLE HYSTEROMETER 30 cm - box of 25**  
The most versatile, economical and eco-friendly single-use hysterometer. Estimating the uterine depth, checking cervical canal permeability and womb orientation. Sterile, disposable medical device, class I. Latex free, PVC free. Shelf life: 5 years. Multilingual manual: GB, FR, IT, ES, PT, DE.

