


## STERILE PAP TEST KITS

Instructions:  
GB, IT, PT

CE 0373




29702

READY TO USE  
PAP TEST KIT

### • 29702 PAP TEST KIT - sterile - box of 100 pcs.

Pap-Test is a painless examination to detect the presence of any anomalous cells on the cervix. It allows an early diagnosis of uterus neck cancer. The functional configuration of the kit allows a ready to go operation.

The kit includes: two latex gloves, one perno speculum (medium size), one brush, one Ayre spatula (wooden), two slides and one box for them. Individually packed in sterile pouch.



29704

### • 29704 PAP TEST KIT - sterile - box of 50 pcs

The kit includes: one pair latex glove, one side-screw speculum (medium size), one GIMA brush ball, one Ayre spatula A size (wooden), two slides. Individually packed in sterile pouch.

### • 29703 PAP TEST KIT - sterile - box of 100 pcs

Same as 29702 but without speculum. Instructions: GB, FR, IT.

29703

## LUBRIGIMA: LUBRICANT GEL

### • 29960 LUBRIGIMA - 250 ml

Facilitates instruments insertion during gyn. examinations. Its special composition and high transparency do not alter the visual image when optical instruments are used. It does not contain spermicides and does not grease or stain.

It has acid Ph and good electrical conductivity. Water soluble.

Well absorbed by the skin, does not irritate the mucosa. LATEX FREE.



29960

## KLY LUBRICATING GEL

### • 29957 KLY LUBRICATING GEL

sachet 5 g - box of 100

### • 29959 KLY LUBRICATING GEL

tube 82 g - 1 tube

### • 29958 KLY LUBRICATING GEL

tube 42 g - sterile and paraben-free - 1 tube  
Hypoallergenic lubricating gel, does not irritate skin, dissolves in water and is easily cleaned. It does not contain oil or oily substances. It has no toxic effects or odour. It is used in endoscopy, gastroscopy, colonoscopy and rectoscopy.



Each 100 g contains  
14 g of glycerine and  
preservatives



29959

Multilingual box and manual (29958-9 only):  
GB, FR, IT, ES, DE, PT, GR, RU, TR, Arabic

## LUBRICATING GEL WITH CHLOREXIDINE

### • 29961 FARMAGEL® C LUBRICATING GEL WITH CHLOREXIDINE - single use squeeze bulb 12.5 g - sterile - box of 25


Farmagel® C is a sterile water-soluble, grease-free, transparent lubricating gel. Suitable for endoscopic procedures which require lubricating action. It reduces friction between probe and mucosa facilitating the entry of diagnostic or therapeutic devices into the body. It does not irritate and is easily removable from skin and clothing.

Ingredients: chlorhexidine (0.05%), hydroxyethyl cellulose, glycerin, purified water, hydrochloric acid. Class IIa medical device.


MADE IN ITALY



29961



## PROCTOSCOPE



### • 25957 PROCTOSCOPE adult - box of 50 pcs.

All proctoscopes are made of PS and packed in single sterile pouch. Disposable product for single use.

mm	A	B	C	D	E
25957	110.5	155	95	17.8	22