

## AYRE SPATULA - BIRCH WOOD

Non sterile 29743

Spatula di ayre sterile in legno  
AYRE SPATULA - WOOD Type A - sterile  
CE 0476 STERILE  
Il utilizzo può causare infezioni crociate Reuse can cause cross-infections

Sterile 29745

GIMA code	AYRE SPATULA - WOOD	Length	Minimum order
29743	Ayre Spatula A 5 boxes of 100 pcs = 500 pcs	18 cm	500 pcs.
29745	Sterile Ayre Spatula A Individually packed in sterile pouch, in box of 100	18 cm	100 pcs.

## AYRE SPATULA - PLASTIC

29748

29749

Used for ecto-cervical and posterior fornix cells sampler in uterus neck cancer screening procedures. 29746

GIMA code	AYRE SPATULA - PLASTIC	Material	Length	Minimum order
29746	Ayre Spatula Packed in plastic bags of 25 pcs., in carton of 400 pcs.	PP	18 cm	400 pcs.
29748	Sterile Ayre Spatula A Individually packed in sterile pouch, in box of 100	PP	18 cm	100 pcs.
29749	Sterile Ayre Spatula B Individually packed in sterile pouch, in box of 100	PP	18.8 cm	100 pcs.

## GIMA BRUSH

Sterile pack 29742

brush Ø=6 mm

CE 0123

MADE IN ITALY

20 cm

CE 29739-29740

- 29739 GIMA BRUSH - conical - bag of 25 pcs.
  - 29740 GIMA BRUSH - conical - box of 400 pcs.
- These single-use endocervical brushes are an efficient method of collecting vulval and endocervical samples for cytological/microbiological examinations. Non sterile. Packed in plastic bags of 25 units, in carton of 400 units.
- 29742 GIMA BRUSH - conical - STERILE - box of 500 pcs
- As code 29740, individually packed in sterile plastic bag. Multilingual manual: GB, FR, IT, ES, PT, DE, FI, GR, HR.

## GIMA BRUSH - STERILE

CE 0476


Special "ball" to avoid cervix puncture and patient pain during the procedure

29736

29735

20 cm

- 29735 GIMA BRUSH B - sterile - box of 100 pcs.
  - 29736 GIMA BRUSH BALL - sterile - box of 100 pcs.
- Each brush is packed in a single sterile pouch and box of 100 pouches.



## COTTON SWABS

29753 Label to write patient data

STERILE SWAB STICK

29754

29756

CE

GIMA code	COTTON SWABS	Shaft Material	Cotton Ø	Shaft length	Minimum order
29750	Cotton Swab Packed in plastic bags of 25 pcs., in carton of 400 pcs. Multilingual instructions: GB, FR, IT, ES, PT	Plastic	1 cm	20 cm	400 pcs.
29753	Sterile Cotton Swab Packed in box of 100 pcs.	Plastic	0.5 cm	15 cm	100 pcs.
29754	Cotton Swab Packed in plastic bags of 100 pcs., in carton of 1000 pcs.	Wood	0.5 cm	15 cm	1,000 pcs.
29756	Cotton Swab Packed in plastic bags of 50 pcs., in carton of 500 pcs.	Wood	1 cm	15 cm	500 pcs.

## GIMA COLLECTOR - STERILE

29737

With Gima collector all it takes is a simple rotating motion to collect a total sampling endocervical and ectocervical cells, including the transformation area, with less trauma. 29737

GIMA code	COLLECTOR	Material	Shaft length	Minimum order
29737	Sterile GIMA collector Packed in a single sterile pouch in box of 100 pouches	Plastic	19.5 cm	100 pcs.