


GISAFE 98% FILTERING 3 PLY SURGICAL MASKS TYPE IIR WITH LOOPS



98% FILTERING 3 PLY MASKS - TYPE IIR with ear loops
 Lightweight and comfortable surgical masks suitable for hospital, dental and outpatient use.

- large range: adult mask in 5 colours, 7 patterns;
- pediatric mask in 1 colour or 1 pattern
- available in box of 50 or flowpack of 10 pcs
- conformable nose band for anatomical fitting
- fiber glass free, filtration efficiency BFE ≥98%
- non woven material, latex free
- 1st ply, outside 28 gsm non woven
- 2nd ply, 25 gsm melt blown filter
- 3rd ply, inside 25 gsm non woven



Multilingual box:
 GB, FR, IT, ES, PT, DE,
 GR, SE, RO, Arabic



Front

Back

20691

TECHNICAL SPECIFICATIONS

BFE%	≥ 98
Differential pressure (Pa/cm ²)	< 60
CfU/g	< 30
Size cm:	17.5x9.5 (adult) - 14x9 (pediatric)
EN 14683:2019 + AC:2019 conform; MDR 2017/745 Class I	

GISAFE 3 PLY SURGICAL MASK - ear loops			
GIMA code	ADULT MASKS - box	Colour / Pattern	Minimum order
20697	3-ply masks - coloured	pink	box of 50
20698	3-ply masks - coloured	white	box of 50
20699	3-ply masks - coloured	light green	box of 50
20700	3-ply masks - coloured	light blue	box of 50
20701	3-ply masks - coloured	black	box of 50
ADULT MASKS in single transparent pouch - box			
20670	3-ply masks - coloured	light blue	box of 50
ADULT MASKS - flowpack			
20694	3-ply masks - coloured	pink	flowpack of 10
20695	3-ply masks - coloured	light blue	flowpack of 10
20696	3-ply masks - coloured	light blue	box of 2000*
20685	3-ply masks - patterned	cartoon	flowpack of 10
20686	3-ply masks - patterned	peace	flowpack of 10
20687	3-ply masks - patterned	stars	flowpack of 10
20688	3-ply masks - patterned	colours	flowpack of 10
20689	3-ply masks - patterned	macarones	flowpack of 10
20690	3-ply masks - patterned	skulls	flowpack of 10
20691	3-ply masks - patterned	military	flowpack of 10
PEDIATRIC MASKS - flowpack			
20680	3-ply masks - coloured	light blue	flowpack of 10
20682	3-ply masks - patterned	cartoon	flowpack of 10

*200 flowpacks of 10 pieces