


DISPOSABLE GROUND PADS - WITH/WITHOUT CABLE



GIMA code	DISPOSABLE GROUND PADS WITH SOLID GEL - need cable 30563	Minimum order
30491	REM GROUND PADS - pediatric	box of 25
30487	REM GROUND PADS - pediatric	box of 100
30492	GROUND PADS - adult	box of 25
30489	GROUND PADS - adult	box of 100
30493	REM GROUND PADS - adult	box of 25
30488	REM GROUND PADS - adult	box of 100
30495	REM NON WOVEN GROUND PAD - ad/ped	box of 100
30563	CABLE FOR PADS - 5 m	box of 1
DISPOSABLE GROUND PADS WITH CABLE - 3m		
30476	PADS WITH CABLE Ø 6.3 mm plug	box of 25
30478	PADS WITH CABLE Valleylab type	box of 25
30479	REM NON WOVEN GROUND PADS with solid gel with cable Valleylab type	box of 25

All plates and cables are individually packed



CONNECTION CABLES FOR DISPOSABLE GROUND PADS

CABLES TO CONNECT GROUND PADS TO ELECTROSURGERY UNITS


The bipolar universal connection cables are supplied with a bipolar clamp for connection to disposable ground pads with a tab (dimensions 22-25 mm).

Various connectors make the cables compatible with numerous electrosurgical devices.

The cable can be fixed to the surgical sheet during the operation thanks to a metal sliding clamp.

Latex free. Individually packed in sealed polyethylene pouch compliance to 93/42/EEC. MD Class I.

30563



COMPATIBLE ELECTROSURGERY UNITS

• 30563 CONNECTION CABLE FOR PADS - 5 m - Universal jack mono plug 6.3 mm

Aesculap: TM, TM401
Alexandra: D205, E300, 500D, 102/201, 2020/E, 202, MICRO180
ALSA: All
 MF100, MF360/A360, 380100, 180
Berchold: Elektrotom 200, 390, 400, 540
Birtcher: 772, 773, 737, 737XL, 774

Bovie: Bovie CSV, URO, BS, 400/400, 400SR, X10, X10L
Burdick: U-7, SU-8
Cameron-Miller: 26-1290
Cobe: TS-403
Concept: 9600, 9700, 9800, 9900
Dolley: COBII 203T, 502, 503T, 505T
Ems Davol: System 2000, 5000

Faset: Serie COBI
Gima: all models
GPS: all models
I Cuatro: Megeri1
Lamidey: Series Subtilec: All
 Series Surgilec: All
Led: all models
Martin: all models

Neomed: Series 3000, 3000A, 3002, 3006; 3010, 3020
Philips: TS-403
Ritter/Mdt: Solide State 400, 400B, 400SR
Sidevan: Seb2
Soxil: BP 2000, X10
Thackray: T300S
Torry: Lem 300 CST

• 30546 CONNECTION CABLE FOR PADS - plug with spacer - 5 m - Universal jack mono plug 6.3 mm


Erbe: ACC, ICC, VIO, T

• 30547 CONNECTION CABLE FOR PADS VALLEYLAB TYPE - 5 m - Universal Valleylab connection

Aspenlabs: Excalibur
Valleylab: REM type devices
Bowa:

Conmed
Erbe ICC International
Erbe ACC International

AUTOCLAVABLE METAL AND RUBBER PLATES



- 30566 METAL PLATE - needs cable 30560
- 30560 CABLE FOR RUBBER STEEL PLATES 30555, 30556, 30566
- 30564 SURGEON NEUTRAL PLATE 16x12 cm - with cable
- 30561 ROD NEUTRAL ELECTRODE with cable for 30625*
- 30562 FLEXIBLE RUBBER PLATE 15x8 cm with cable 4.5 m
- 30565 FLEXIBLE RUBBER PLATE 25x15 cm with cable 4.5 m

*Available on request for other electrocautery

- 30555** FLEXI STEEL RUBBER PLATE 18x21 cm - needs cable 30560
- 30556*** FLEXI STEEL REM RUBBER PLATE - needs cable 30560
- 30560 CABLE FOR RUBBER STEEL PLATES 30555, 30556, 30566
- 30574 REM RUBBER PLATE 20x15 cm - needs cable 30575
- 30575 CABLE FOR RUBBER PLATES 30573 and 30574

**Compatible with all Gima diatermo

***Compatible with all Gima diatermo except 30516, 30625