



## BIPOLAR CABLES FOR MAJOR BRANDS



Quality highly flexible silicone bipolar cables 3 m long (10").  
Autoclavable at 134°C for minimum 20 times.  
Recommended EtO sterilization for max extended life.

Suitable for most generators, with connector for European or US  
forceps. Multilingual box and manual: GB, FR, IT, ES, PT, DE, SE PL,  
RO,GR, HR, Arabic

|                                                                                                                                                                        |                                                                                                           |                                                                                      |                                                                                                                                                                              |                                                                                      |                                                                                                           |                                                                                      |                                                                                                                                              |                                                                                        |                                                                                                           |                                                                                        |                                                                                                                                        |  |                                                                                                             |  |                                                                                                                          |  |                                                                                                             |  |                                                                                                                          |  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------------------------------------------|--|--------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------------------------------------------|--|--------------------------------------------------------------------------------------------------------------------------|--|
| <br><b>30642</b><br>• <b>30642</b><br><b>BIPOLAR</b><br><b>CABLE</b><br><b>ADAPTOR</b> | <br>Connector to forceps |                                                                                      | <br>Connector to electrosurgical unit                                                       |                                                                                      | <br>Connector to forceps |                                                                                      | <br>Connector to electrosurgical unit<br>Pin distance 29 mm |                                                                                        | <br>Connector to forceps |                                                                                        | <br>Connector to electrosurgical unit<br><b>30642</b> |  | <br>Connector to forceps |  | <br>Connector to electrosurgical unit |  | <br>Connector to forceps |  | <br>Connector to electrosurgical unit |  |
|                                                                                                                                                                        | GIMA Code                                                                                                 | <b>30632</b>                                                                         | <b>30618</b>                                                                                                                                                                 | <b>30640</b>                                                                         | <b>30637</b>                                                                                              | <b>30643</b>                                                                         | <b>30618+ 30642</b>                                                                                                                          | <b>30617</b>                                                                           | <b>30644</b>                                                                                              | <b>30639</b>                                                                           | <b>30638</b>                                                                                                                           |  |                                                                                                             |  |                                                                                                                          |  |                                                                                                             |  |                                                                                                                          |  |
| Connector to forceps                                                                                                                                                   |  EU                      |  US |  EU                                                                                         |  US |  EU                      |  US |  EU                                                       |  US |  EU                    |  EU |                                                                                                                                        |  |                                                                                                             |  |                                                                                                                          |  |                                                                                                             |  |                                                                                                                          |  |
| Compatibility brand/model                                                                                                                                              | GIMA: MB 240<br>MB 380<br>VALLEYLAB,<br>CONMED,<br>CODMAN,<br>OLYMPUS                                     |                                                                                      | GIMA: MB 120 F, MB 200 F,<br>MB 200 D, MB 200T,<br>MB 300 D, MB 400 D,<br>MB 300T, MB 400T<br>LED: SURTRON 120, 200, 200T,<br>300, 300T, 400, 400T<br>VALLEYLAB, EMC, DOLLEY |                                                                                      | GIMA: MB 122, MB 160,<br>MB 200, MB 202                                                                   |                                                                                      | GIMA: MB 80 D, MB 120 D,<br>MB 160 D<br>LED: SURTRON 80, 120,<br>160                                                                         |                                                                                        | MARTIN<br>BERCHTOLD<br>AESCULAP<br>WOLF                                                                   |                                                                                        | ERBE<br>SELECT<br>DOWNS<br>SIEMENS                                                                                                     |  |                                                                                                             |  |                                                                                                                          |  |                                                                                                             |  |                                                                                                                          |  |



## DISPOSABLE BIPOLAR FORCEPS - VALLEYLAB CABLE

- 30300 DISPOSABLE JEWELER FORCEPS - 11.5 cm - tip 0.5 mm
- 30302 DISPOSABLE STRAIGHT FORCEPS - 20 cm - tip 1.0 mm
- 30304 DISPOSABLE STRAIGHT FORCEPS - 15 cm - angled tip 1.0 mm
- 30306 DISPOSABLE BAYONET FORCEPS - 20 cm - tip 1.0 mm

Disposable bipolar forceps are supplied with 3 m cable.

Valleylab plug (same connector and compatibility as 30640).

Sterilization by Ethylene Oxide (EtO).



## BIPOLAR SCISSORS


- 30310 BIPOLAR SCISSORS - 18 cm (7") - curved
- Needs bipolar cable 30315
- 30315 SILICONE BIPOLAR CABLE - 3 m (10')

with Ethicon type connector

For bipolar scissors. Connector 2.5 mm

Compatibility: same models as 30632, 30640, 30617\*, 30643\* (see above).

\*Need adaptor 30642




## BIPOLAR ARTERY SEALER

- 30320 BIPOLAR ARTERY SEALER - 26 cm (10.5")
- Needs bipolar cable 30325
- 30325 SILICONE BIPOLAR ARTERY SEALER CABLE - 3 m (10')

Connector 4 mm, Banana plug.

Compatibility: same models as 30632, 30640, 30617\*, 30643\* (see above).

\*Need adaptor 30642



## FOOTSWITCHES AND BOARDS

| GIMA code                                                                             | FOOTSWITCHES FOR 50, 80, 106, MB 80D, 120D, 120F, 122, 132, 160, 160D, 200, 200D, 200F, 200 T, 202 |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| 30577                                                                                 | MB waterproof footswitch - 2m cable                                                                |
| 30578                                                                                 | MB waterproof footswitch - 3m cable                                                                |
| ANTIEXPLOSION FOOTSWITCH FOR MB 120D, 120F, 122, 132, 160, 160D, 200, 200D, 200F, 202 |                                                                                                    |
| 30602                                                                                 | Double antiexplosion, water proof MB footswitch                                                    |
| ANTIEXPLOSION FOOTSWITCHES FOR MB 240, 300D, 300T, 380, 400, 400D, 400T               |                                                                                                    |
| 30603                                                                                 | Double antiexplosion, water proof footswitch - MB240-380                                           |
| 30604                                                                                 | Double antiexplosion, water proof footswitch - MB250-400-300D-400D                                 |
| 30608                                                                                 | Double antiexplosion, water proof footswitch - MB300-400 TOUCH                                     |
| MAIN BOARD                                                                            |                                                                                                    |
| 30557                                                                                 | Main board - specify model and serial number of diatermo with order                                |
| 30558                                                                                 | Output board for MB 240, 250, 380, 400 - specify as above                                          |



It allows selection from pedal of CUT and COAG functions  
Antiexplosion and waterproof.