# CHANGE-A-TIP™ CAUTERIES WITH HIGH TEMPERATURES

Improved battery-operated cautery with a safety-enhanced design and a recessed button to reduce the potential


30590 for inadvertent activation. (Fine tip, not sterile, included)

It also allows the practitioner to "snapped off" the end of the cautery for safe and effective battery removal.

- interchangeable sterile tips (non sterile tips for 30584 and 30590)
- replaceable batteries
- reusable handles

- no sterilization necessary.





**30584** High temp cautery kit (1 handle + 4 non sterile tips)


## REUSABLE HIGH TEMPERATURE CAUTERY WITH REPLACEABLE TIP

### **CAUTERIES - 1,200°C (2,200°F)**

- interchangeable tips with nickel-chrome wire MADE IN ITALY
- replaceable batteries
- reusable handles

Wide application in gynaecology, orthopedy, ear/ nose/throat, dermatology and urology.

GIMA	INTERCHANGEABLE TIP CAUTERIES
code	MINIMUM ORDER: 10 TIPS
30598	Cautery handle "AA" (tip not included)
30599	High temp long fine tip - sterile - box of 10
30600	High temp large tip - sterile - box of 10
30601	High temp fine tip - sterile - box of 10



Lenath 18 cm


# DISPOSABLE HIGH/MEDIUM/LOW TEMPERATURE CAUTERIES

30701

30700

- 30700 ELECTROCAUTERY 1,200°C fine tip box of 10 Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels. With a fine tip required in most of the surgeries.
- 30701 ELECTROCAUTERY 1,200°C large tip box of 10
- Same as 30700, but with a large tip.
- 30702 LONG SHAFT ELECTRO CAUTERY 1,200°C fine tip box of 10 Same as 30700 but with the tip fixed on a 12 cm long shaft. Precise cauterization even in difficult access sites.
- 30703 EMERGENCY ELECTROCAUTERY 800°C thick tip box of 10 Same as 30700 but with "thick tip" for the perforation of nails in case of haematoma.
- 30704 ELECTROCAUTERY 600°C fine tip box of 10 Same as 30700 but with lower temperature 600°C. Expecially designed for ophthalmology. C 6 2797 MADE IN ITALY

					1
Applications	30700	30701	30702	30703	30704
Cardiovascular, thoracic surgery	Х				
Proctology	Х	Х	Х		
General surgery			Х		
Gynaecology	Х				
Surgical gynaecology			Х		
Neurosurgery	Х	х	Х		
Orthopedic	Х	Х			
Urology	Х				
Obstetrics	Х				
Emergency ward				Х	
Otolaryngology	Х	Х	Х		
Plastic and reconstructive	Х	Х			
Ophtalmology	Х				Х
Dermatology	Х				
Dental surgery	Х				



Features	HTC (High Temperature Cautery)			MTC (Med. Temp. Cautery)	LTC (Low. Temp. Cautery)
	30700	30701	30702	30703	30704
Tip shape	Fine	Large	Fine	Thick	Fine
Tip length	28mm	28mm	125mm	28mm	28mm
Temperature °C/°F	1200/2200	1200/2200	1200/2200	800/1400	600/1100
Weight	65g	65g	78g	65g	45g

Long life

(60 min.)

internal battery