

# VALUE&SAFE FFP2 FILTERING MASKS

## VALUE&SAFE FFP2 FILTERING MASKS WITH EAR LOOPS

Five-layer FFP2 mask to protect yourself and others effectively. Meeting FFP2 requirements of EN 149:2001+A1:2009 in accordance with the PPE regulation EU 2016/425.

- ergonomic design fit face closely
- non irritating, soft and skin-friendly inside non woven fabric
- adjustable nasal clip
- stretchable ear loops, painless after long-time wearing
- foldable in half for easy gathering
- each mask is individually packed in a transparent pouch
- latex free (adult only)

**FFP2 FILTERING EFFICIENCY >94%**

CE 0370



20732, 20733

### FILTERED PARTICLES



VALUE&SAFE FFP2 FILTERING MASKS				
GIMA code	Model	Box and instructions languages	Pattern or colour	Minimum order
<b>FFP2 WITHOUT VALVE</b>				
20738	FFP2 - adult-small	GB, FR, IT, ES, DE	rainbow	box of 20
20739	FFP2 - adult-small	GB, FR, IT, ES, DE	dinosaurs	box of 20
20743	FFP2 - adult-small	GB, FR, IT, ES, DE	white	box of 20
<b>FFP2 WITHOUT VALVE</b>				
20731	FFP2 - adult	GB, FR, IT, ES, DE	light grey	box of 20
20732	FFP2 - adult	GB, FR, ES, DE, PT	white	box of 20
20733	FFP2 - adult	IT, GR, PL, CZ, RO	white	box of 20
20734	FFP2 - adult	IT, GR, DE, SE, PL	blue	box of 20
20735	FFP2 - adult	GB, FR, IT, ES, PT	black	box of 20
20736	FFP2 - adult	GB, IT, PT, RO, CZ	burgundy	box of 20
<b>FFP2 WITH VALVE</b>				
20737	FFP2 - adult	GB, FR, IT, ES, PT	white	box of 20
<b>TECHNICAL SPECIFICATIONS</b>				
Inhalation:	30 l/min% ≤ 0.7	Size:	Adult-small	Adult
Inhalation:	95 l/min% ≤ 2.4	unfolded:	18x13.5 cm	21x15.5 cm
Exhalation:	160 l/min% ≤ 3.0	folded:	9x13.5 cm	10.5x15.5 cm
Penetration of filter material:	≤ 6%			



20738

20739

20743

CE 0598



20731

20734

20736

20735

CE 0370



CE 0598

FFP2 with valve

20737

Each mask is individually packed

# EFFECTIVE FFP2 PARTICLE FILTERING MASKS

• 20740 EFFECTIVE FFP2 FILTERING MASK GB, FR, ES, PT, DE - box of 20 (4 bags of 5)

• 20741 EFFECTIVE FFP2 FILTERING MASK IT, SE, GR, RO, ARABIC - box of 20 (4 bags of 5)


Six-layers filtering FFP2 mask to protect yourself and others effectively.

Meeting FFP2 requirements of EN 149:2001+A1:2009 in accordance with the PPE regulation EU 2016/425.

- ergonomic design fit face closely
- non irritating, and soft inside non woven fabric
- adjustable nasal clip
- flexible and stretchable ear loop, painless after long-time wearing
- foldable in half for easy gathering

Includes 20 hooks for elastic adjustment behind head.

**6 LAYERS**



**FFP2 >94%\***

\*Filtering efficiency



Fit seamlessly Anti-fog

Hook for elastic adjustment

Foldable in half for easy gathering

Elastic ear loop



20740

CE 0598  
FDA



TECHNICAL SPECIFICATIONS	
Inhalation: 30 l/min% ≤ 0.7	Penetration of filter material: ≤ 6%
Inhalation: 95 l/min% ≤ 2.4	Size unfolded: 21.4x16 cm
Exhalation: 160 l/min% ≤ 3.0	folded: 10.7x16 cm