


## E-SCOPE®: RIESTER OTO-OPHTHALMOSCOPES

- 31546 E-SCOPE® OTOSCOPE - HALOGEN 2.7V - white in pouch
  - 31548 E-SCOPE® OTOSCOPE - HALOGEN 2.7V - black in pouch
  - 31549 E-SCOPE® OTOSCOPE - XENON 2.5V - white in pouch
  - 31550 E-SCOPE® OTOSCOPE - XENON 2.5V - black in pouch
  - 31552 E-SCOPE® F.O. OTOSCOPE - LED 3.7V - white in rigid case
  - 31553 E-SCOPE® F.O. OTOSCOPE - LED 3.7V - black in rigid case
- Halogen, direct Xenon or F.O. LED illumination otoscopes.  
Sealing precision optics pivoting to both sides with triple

3X magnifying pivoting lens



magnification for superior view

Suitable for performance of pneumatic tests (supplied without connector and insufflator)

- 31556 E-SCOPE® OPHTHALMOSCOPE - XENON 2.5 V - black - in rigid case
- 31557 E-SCOPE® OPHTHALMOSCOPE - F.O. - LED 3.7 V - white - in rigid case
- 31559 E-SCOPE® OPHTHALMOSCOPE - F.O. - LED 3.7 V - black - in rigid case

E-scope ophthalmoscopes are available with Xenon or LED illumination.

- with 2.5V Xenon or innovative 3.7V LED illumination 31550 for higher-contrast and more efficient diagnosis
  - dioptr disc with 18 corrective lenses ( $\pm 1, 2, 3, 4, 6, 8, 10, 15, 20$ )
  - easy-to-use aperture wheel with six apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
  - high performance optics with aspherical condenser
  - spectacle protection
  - parallel optical path
- MADE IN GERMANY



- 31558 E-SCOPE® DIAGNOSTIC KIT - white
  - 31564 E-SCOPE® DIAGNOSTIC KIT - black
- Oto-ophthalmoscope sets, LED 3.7V, in rigid case.
- 31560 E-SCOPE® OTOSCOPE - XENON 2.5 V BULB
  - 31561 E-SCOPE® OTOSCOPE - LED 3.7 V BULB
  - 31562 E-SCOPE® OPHTHALMOSCOPE - XENON 2.5 V BULB
  - 31563 E-SCOPE® OPHTHALMOSCOPE - LED 3.7 V BULB
  - 31817 RIESTER® EAR SPECULUM Ø 2.5 mm - autoclav. - box of 10
  - 31819 RIESTER® EAR SPECULUM Ø 4 mm - autoclav. - box of 10
  - 31566 E-SCOPE® EAR SPECULUM Ø 2.5 mm - disposable - box of 100
  - 31567 E-SCOPE® EAR SPECULUM Ø 4 mm - disposable - box of 100

## LUXAMED AURIS LED OTOSCOPES



- 31570 LUXAMED AURIS LED 2.5V OTOSCOPE - black
  - 31571 LUXAMED AURIS LED 2.5V OTOSCOPE - blue
  - 31572 LUXAMED AURIS LED 3.7V USB OTOSCOPE - black
- Innovative otoscopes that combine best handcraft production with selected materials with the most modern technology. Two-components handle (aluminium/fibre-glass reinforced) with electronic ON/OFF switch and dimmer.

Features:

- LED-Ring illumination with 6 LEDs (concentric), no fibre optic. The functionality is not affected (patent DE202013104278U1) if a single LED fails.

**LUXAMED®**

- automatic switchoff after 3 minutes
- approx. 100,000 hours life expectancy of LEDs
- approx. 10,000 Lux intensity of illumination
- approx. 4,000 K colour temperature
- pivoting 3-fold magnifying glass with locking function
- insufflation port for pneumatic test.

Supplied with 2xAAA alkaline batteries (31570/1) or lithium-ion rechargeable battery 3.7V (31572) 20 disposable ear funnels (10x2.5 mm and 10x4.0 mm) and manual (GB, IT, DE). See page 266 for compatible ear speculum.

MADE IN GERMANY



Recharger with USB plug - 3 adapters (EU, UK, US)