

RIESTER ELITEVUE



31863

Riester
MADE IN GERMANY

31860

Supplied in a case with specula

- 31860 RIESTER ELITEVUE LED OTOSCOPE 3.5 V - 2200-203
 - 31863 RIESTER ELITEVUE LED OTO-OPHTHALMO 3.5 V - 2210-203
- Both with C-handle for rechargeable Li-Ion battery**
EliteVue offers twice the field of view and 66% greater magnification (at 15 mm distance) compared to standard otoscopes. Improved view of auditory canal and tympanic membrane guarantees more efficient diagnosis. Suitable for throat illumination. 4,000 K LED illumination (CRI: 92). 5X magnification.
- 31840 RI-CHARGER L 230 V - 10705
 - 31841 RI-ACCU L 3.5 V Li-Ion BATTERY for "C" size - 10691

HEINE K180 OTO-OPHTHALMO



5 YEARS WARRANTY



MADE IN GERMANY

31740

- 31740 HEINE K180 - F.O. OTO-OPHTHALMO SET - 2.5 V - A-279.10.118
- Otoscope and ophthalmoscope K180 with BETA battery 2.5 V handle and a set of single use specula. Supplied in hard case.
- 31782 HEINE 084 OPHTHALMO BULB 2.5 V - spare - X-001.88.084
 - 31772 HEINE 077 OTOSCOPE BULB - spare - X-001.88.077
 - 31741 HEINE K180 OPHTHALMOSCOPE - 2.5 V C-182.10.118
- 28 lenses from -35D to +40D supplied in pouch



HEINE BETA 200 F.O. OTOSCOPES



MADE IN GERMANY

31725 Large, scratch-resistant glass viewing window with 3x magnification

34449

Handle with charge status display

31728+34374

Distal fiber optic (F.O.) illumination with XHL xenon halogen technology
Airtight window seal for reliable test of tympanic mobility

- 31725 HEINE BETA 200 F.O. OTOSCOPE - 2.5V - B-141.10.118
- F.O. Otoscope with BETA battery 2.5V handle with rheostat, a set of single use specula, a set of reusable specula and a spare bulb. Supplied in hard case. Solid all metal construction, a quality instrument with fibre optic illumination.
- 31772 HEINE 077 OTOSCOPE BULB - spare - X-001.88.077
 - 31784 HEINE 078 OTOSCOPE BULB - X-002.88.078

- 31728 HEINE BETA 400 LED F.O. OTOSCOPE - B-143.24.420
- Exclusive 4.2X magnification. LED in HQ with continuous brightness control. All metal housing and scratch resistant glass lenses. Colour temperature: 3,500 K. Colour rendering index (CRI) >97. Exceptional depth of field, provides sharp spatial image. LED otoscope with rechargeable handle supplied with 4 reusable speculum, 5 of both All Spec disposable speculum 2.5 and 4 mm Ø in a case with NT4 table charger (34374).

HEINE BETA 200 AND 200S OPHTHALMOSCOPES



31730

6 standard diaphragms: slit, fixation star with polar co-ordinates, cobalt blue filter, large spot, small spot, hemi spot

6 diaphragms with green filter

+ in 1 D steps 1 - 10 15 20 40 D
- in 1 D steps 1 - 10 15 20 25 35 D

- 31730 HEINE BETA 200 HALOGEN OPHTHALMOSCOPE - 2.5 V - C-144.10.118
- The unique aspherical optic system eliminates corneal and iris reflexes. The illuminated area of the retina remains fully visible even in case of small pupil (bright light). BETA 200 has 6 diaphragms with separately "green filter". It is recommended for small or expanded pupils; suitable for slit or fluorescent ophthalmoscopy as well as for fixation test. Range -35D to +40D. Supplied in hard case with replacement bulb.
- 31731 HEINE BETA 200 F.O. OPHTHALMOSCOPE - 3.5 V Li-ion - C-144.27.388
- Same as 31730 but with USB rechargeable handle for mains socket and Li-ion batteries.
- 31783 HEINE 070 OPHTHALMOSCOPE BULB 3,5 V - X-002.88.070



31730

- 31735 HEINE BETA 200S OPHTHALMOSCOPE - 2.5 V - C-261.10.118
- 31781 HEINE 069 OPHTHALMO BULB - 2.5 V - spare - X-001.88.069
- 31737 HEINE BETA 200S LED OPHTHALMO - 3.5 V with USB rechargeable handle - C-261.28.388

Superior aspherical optics and 74 single-diopter steps - optical system: corneal and iris reflexes are eliminated. A full view is guaranteed even with small pupils.

- XHL (Xenon Halogen Technology) or LED: bright, white light.
- recessed, multi-coated viewing window: avoids stray light.
- optics mounted on precision metal chassis.
- ergonomic shape: instrument fits the orbita comfortably in any position.
- 7 apertures with separate red-free filter:
- lens range: perfect focus even with high refractive errors.
- exact colour rendering both with XHL and LED illumination

MADE IN GERMANY



Cobalt blue filter, fixation star with polar co-ordinates, large spot, small spot, pinhole, slit, hemisport

With red-free filter: Cobalt blue filter, fixation star with polar coordinates, large spot, small spot, pinhole, slit, hemisport

74 single diopters steps
+ in 1 D steps: +1 through +38 D
- in 1 D steps: -1 through -36 D

