

DISPOSABLE EAR-SPECULUM FOR "TRADITIONAL" OTOSCOPES

GIMA code	DISPOSABLE EAR - SPECULUM FOR OTOSCOPES OF DIFFERENT BRANDS	Colour	To be used with	Minimum order
31315	Disposable ear-speculum Ø 4 mm	Black	Riester (some models)	box of 200 pcs.
31504	Disposable ear-speculum Ø 2.5 mm	Black	Heine Beta 100, K100 and Eurolight,	box of 1000 pcs. (20 pack of 50)
31505	Disposable ear-speculum Ø 4 mm	Black	Kawe C10, C30, Combilight C10	box of 1000 pcs. (20 pack of 50)
31506	Disposable ear-speculum Ø 2.5 mm	White	Riester: Ri-star	dispenser of 250 pcs.
31507	Disposable ear-speculum Ø 4 mm	White		dispenser of 250 pcs.
31508	Disposable ear-speculum Ø 2.5 mm	Black	Welch Allyn: Pocket-pro; Pocket scope	dispenser of 250 pcs.
31509	Disposable ear-speculum Ø 4 mm	Black	Series 200, 250	dispenser of 250 pcs.

Suitable for RIESTER  31315	Suitable for HEINE  31504-31505	Suitable for RIESTER Multilingual box: GB, FR, IT, ES, PT, DE, GR, DK, Arabic.  31506-31507	Suitable for WELCH-ALLYN Multilingual box: GB, FR, IT, ES, PT, DE, GR, DK, Arabic.  31508-31509
--	--	--	---

WELCH-ALLYN KLEENSPEC DISPOSABLE EAR SPECULUM

GIMA code	WELCH-ALLYN KLEENSPEC Disposable ear-speculum	Colour	To be used with	Minimum order	
72570	KleenSpec ear-speculum Ø 3 mm	Black	Welch-Allyn 20100, 20200 20202, 21600, 21700	box of 100 (4 packs of 25 pcs.)	 72570
32105	KleenSpec ear-speculum Ø 2.75 mm	Black	Welch-Allyn, Macroview	box of 850	
32106	KleenSpec ear-speculum Ø 4.25 mm	Black	Riester Ri-scope L2, L3	box of 850	

VARTA ALKALINE, LITHIUM AND HEARING AID BATTERIES

								
32200	32201	32203	32204	32205	32216	32217	32218	32219

FOR GENERAL USE APPLICATIONS **MADE IN GERMANY** **FOR ACOUSTIC**

GIMA code	VARTA ALKALINE LONG LIFE POWER BATTERIES	Type	Blister	GIMA code	VARTA "HEARING AID BATTERIES" - mercury free	Blister
32200	Varta ministilo - alkaline	AAA	4 pcs.	32216	Varta 675 - acoustic	6 pcs.
32201	Varta stilo - alkaline	AA	4 pcs.	32217	Varta 13 - acoustic	6 pcs.
32203	Varta half torch - alkaline	C	2 pcs.	32218	Varta 312 - acoustic	6 pcs.
32204	Varta torch - alkaline	D	2 pcs.	32219	Varta 10 - acoustic	6 pcs.
32205	Varta 9V - alkaline*	9V	1 pc.			

*Made in China



		
32210	32211	32212

FOR ELECTRONIC

GIMA code	VARTA INDUSTRIAL ALKALINE BATTERIES	Type	Box	GIMA code	VARTA LITHIUM BATTERIES - BUTTON	Blister
32226	Varta Industrial ministilo - alkaline	AAA	10 pcs.	32210	Varta 2016 - lithium button	1 pc.
32228	Varta Industrial stilo - alkaline	AA	10 pcs.	32211	Varta 2025 - lithium button	1 pc.
				32212	Varta 2032 - lithium button	1 pc.

RECHARGEABLE BATTERIES AND CHARGERS

GIMA code	VARTA POWER PLAY RECHARGEABLE BATTERIES	Type	Blister
32235	Varta ministilo - 1,000 mAh	AAA	2 pcs.
32236	Varta stilo - 2,600 mAh	AA	4 pcs.

GIMA code	CHARGERS FOR ANY BRAND**	Power supply	Sizes mm	Weight g	To be used with:	Recharging time	12V+USB adaptor
32241	Universal charger	100-240 V	172x109x53	403	AAA, AA, C, D, 9C	5 h	no
32238	LCD charger	110-240 V	72x116x32	170	AAA, AA	2 h	yes

**Supplied without batteries

			
32235	32236	32238	32241