

OTOCLEAR® PORTABLE ELECTRONIC IRRIGATOR

• 25855 OTOCLEAR® PORTABLE ELECTRONIC IRRIGATOR

The OtoClear® Portable Electronic Irrigator is a versatile solution for automated ear irrigation. Designed for safe and effective cerumen removal when used with patented OtoClear Tips. The Tip's unique, flared design prevents over-insertion while prohibiting any backsplash. The three-stream water flow creates a circular clearing method around the ear canal wall, directing water away from the tympanic membrane.

The result is a safe, simple and comfortable approach to ear wax removal. Recommended for high volume facilities.

Features:

- adjustable flow settings
- effective, pulsed irrigation
- electronic pressure-control system with three settings: low, medium, high
- 278 ml self-contained water reservoir
- rechargeable battery

Supplied with: ear basin, 5 OtoClear tips and 220V charging adapter.



• 25861 OTOCLEAR® EAR IRRIGATION TIPS - single use - box of 40 - spare

OtoClear® Ear Irrigation Tip is a soft, flared, flexible, 3-streamed-tip that improves patient safety, prevents the risk of over-insertion and eliminates the mess of cerumen removal during irrigation. Flared design prevents backplash and over-insertion. It directs water to the ear canal wall, with three divergent streams instead of directly at the tympanic membrane. Compatible with a variety of delivery devices.

• 25856 ADAPTOR WAND - spare

• 25857 WATER RESERVOIR 278 ml - spare

• 25862 EAR IRRIGATION BASIN

- box of 3 - spare 25855

MADE IN THE USA



GLASS-METAL AND ALL-METAL EAR SYRINGES



GLASS-METAL JANET RECORD SYRINGES

High quality, reliable and accurate glass-metal syringes. Janet Record come complete with 2 cannulas, one conical, one oval. Individually packed in cardboard box.

GIMA code	GLASS-METAL SYRINGES COMPLETE WITH 2 CANNULAS
25801	Janet Record 50 cc - 13 x Ø 3.5 cm
25802	Janet Record 100 cc - 14 x Ø 4.7 cm
25803	Janet Record 200 cc
25804	Spare cannula for Janet Record (conical shape)
25805	Spare cannula for Janet Record (olive shape)
25861	Ear irrigation tip - single use - box of 40




ALL-METAL SCHIMMELBUSCH SYRINGES

High quality, reliable and accurate syringes. Schimmelbusch come complete with 3 cannulas and ear bowl. Individually packed in cardboard box.

GIMA code	ALL-METAL SYRINGES COMPLETE WITH 3 CANNULAS AND EAR BOWL
25806	Schimmelbusch 50 cc
25807	Schimmelbusch 100 cc
25808	Schimmelbusch 200 cc
25809	Spare cannula for Schimmelbusch (conical shape)
25812	Spare cannula for Schimmelbusch (olive shape) small
25813	Spare cannula for Schimmelbusch (olive shape) large
25861	Ear irrigation tip - single use - box of 40

HIGH CAPACITY POLYPROPYLENE SYRINGES



HIGH CAPACITY POLYPROPYLENE SYRINGES - sterile

Sterile, single use, latex free syringes for ear and other purposes. Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.

Medical Device class IIa.

GIMA code	HIGH CAPACITY, POLYPROPYLENE SYRINGES STERILE - CATHETER CONE	Minimum order
23820	Syringe 60 ml	box of 20
23821	Syringe 100 ml	box of 25