

## OTOCLEAR® EAR IRRIGATION SYSTEM

Kit includes:  
1 spray wash bottle  
20 otoclear tips  
1 ear basin



25858



- **25858 OTOCLEAR EAR SPRAY WASH**  
- manual, intermittent lavage  
- delivers greater volume compared to syringes or other manual devices  
- efficient design allows for two handed procedure

MADE IN THE USA


- **25861 OTOCLEAR EAR IRRIGATION TIPS**  
single use - box of 40 - spare  
Flared design protects against over-insertion. The tip directs three streams of water towards the walls of the ear canal, away from the tympanic membrane, creating a turbulent lavage that effectively dislodges the cerumen. The exit portal controls backsplash and directs drainage into a basin.
- **25862 EAR IRRIGATION BASINS** - box of 3 - spare



Thermometer for water temperature



Improved patient comfort thanks to the soft "gentle touch" design



Patented 3-stream action

## LOOP-SPOON EAR CURETTES

- **31470 GIMA EAR CURETTES** - box of 10 pcs  
Double function currettes: one side spoon, one side loop.  
Each curette is individually packed in a single pouch.

31470



## BIONIX EAR CURETTES



**FLEXLOOP:**  
For normal curettage  
Most flexible of all  
Least traumatic  
Rounded tips for comfort  
4 mm wide tip



**MICROLOOP:**  
For compact cerumen  
More rigid construction  
Tip can be angled for added strength  
Rounded tip for comfort  
3 mm wide tip



**INFANTScoop:**  
For the smallest ears  
For when a small scoop works better  
Can be used through an otoscope  
Flexible 2 mm wide tip



**CERASPOON:**  
For viscous cerumen  
For extensive cerumen buildup  
For when a large scoop works better  
4 mm wide tip

The complete cerumen management system

- **31462 FLEXLOOP** - white - box of 50 pcs
- **31464 INFANTScoop** - light blue - box of 50 pcs
- **31460 MICROLOOP** - green - box of 50 pcs
- **31467 CERASPOON** - yellow - box of 50 pcs


Our four styles of currettes give you the most effective, economical and safe way to remove cerumen - whatever the size of ear or consistency of cerumen.

**Patient comfort:** Metal currettes are cold, hard and sharp. BIONIX ear currettes with rounded plastic tips reduce trauma and risk of laceration. Latex free.

**Effectiveness:** Ear currettes can duplicate the effectiveness of metal currettes in almost any situation and their polypropylene composition makes them virtually unbreakable.

**Disposability:** Ear currettes don't require expensive and time-consuming sterilization and eliminate the risk of cross-contamination.

MADE IN THE USA



## LIGHTED EAR CURETTES

- **31450 KIT MIXED LIGHTED CURETTES**  
Kit is composed by 50 lighted currettes  
(8 FlexLoop and 7 of each of below other 6 tip styles) + LED illuminator + magnification lens

MicroLoop®  
3 mm diameter loopFlexLoop®  
4 mm diameter loopInfantScoop®  
2 mm diameter tipAngleLoop®  
Angled 4 mm diameter loopWaveCurette™  
new wave like 2.5 mm tipVersaLoop®  
3 mm wide tipCeraSpoon®  
4 mm wide tip

MADE IN THE USA

- **31461 KIT MICROLOOP**
- **31463 KIT FLEXLOOP**
- **31465 KIT INFANTScoop**
- **31466 KIT ANGLELOOP**
- **31468 KIT VERSALoop**

All kits are composed by 50 lighted currettes + LED light source + magnification lens

- **31458 LIGHTED FORCEPS** - box of 10 + magnification lens + LED illuminator

- **31459 LIGHTED ARTICULATING CURETTES** - box of 25 + magnification lens + LED light source



Lighted forceps for foreign body removal with LED illuminator



Lighted ear curette 3 mm tip with LED illuminator

Ear curette tip has a 90° range of motion. It articulates to assist in removal of cerumen from the ear canal

