


## GERMAN-MADE DIAPASONS - TUNING FORKS



### ENT TUNING FORKS:

All made in Germany, with age-long craftsmanship experience, in special alloy that ensures maximum sound precision.

GIMA code	ENT TUNING FORKS
31230	Tuning fork 128 Hz* - with weight
31231	Tuning fork 256 Hz* - with weight
31232	Tuning fork 512 Hz
31233	Tuning fork 1024 Hz
31234	Tuning fork 2048 Hz
31235	Tuning fork 4096 Hz
31240	Set of 5 tuning forks in a stylish lined case (contains 128 Hz*, 256 Hz*, 512 Hz, 1024 Hz, 2048 Hz). <i>Tuning forks with asterisk * are supplied with weight</i>

### ENT/NEUROLOGY TUNING FORKS


All made of a special alloy and supplied with base. All tuning forks excluding 31246 have multilingual manual and box: GB, FR, IT, ES, PT, GR. Additional DE, HU, RO, BG only for 31247. GB, IT for 31246. **MADE IN GERMANY**

GIMA code	ENT/NEUROLOGY TUNING FORKS
31245	Tuning fork a1-440 Hz - ENT - 90 g
31246	Tuning fork - LUCAE - 85 g adjustable from c128 to h240 with foot
31247	Tuning fork C64, c128 - RYDEL SEIFFER detachable slider either C64 or c128 - 130 g

## ALUMINIUM ALLOY DIAPASONS - TUNING FORKS

- 31225 ALUMINIUM TUNING FORK 128 Hz
  - 31226 ALUMINIUM TUNING FORK 128 Hz - with foot
- Made of finest aluminium alloy and precisely shaped to meet highest requisites (standards) in medical field.

Resting pad for forehead



## NYSTAGMUS SPECTACLES - DR. FRENZEL

**MADE IN GERMANY**

- 31204 NYSTAGMUS SPECTACLES - DR. FRENZEL - black
- Easy to use in the search of spontaneous or provoked nystagmus. Fixed, biconvex lens + 18.0 diopters. Nystagmus spectacles should be connected only with < 4.2 V transformer as per safety norms 93/42/EEC. See code 31203 or 31214

- 31214 BATTERY HANDLE for 31204
- Offers mobility during examination. Needs 3x1.5 V CLR14 batteries (not included). Manual GB, IT, SE.

- 31203 PLUG-IN TRANSFORMER for 31204 - 100-240 V, 47-63 Hz
  - 31206 BATTERY NYSTAGMUS SPECTACLES - DR. FRENZEL - black
- This model allows wireless working and is provided with power through the adjustable battery placed above the spectacle. Requires three 1.5 Volt batteries, size C LR14/AM-2 (baby cells). Illumination can be adjusted through a turn of the switch.

- 31208 NYSTAGMUS BULB - spare for 31204, 31206
  - 31207 LED BATTERY NYSTAGMUS SPECTACLES - white
  - 31209 LED POWER SUPPLY NYSTAGMUS SPECTACLES - white
- Same as 31206 with LED and inner function that makes easier watching the rashes within a particular time period (30, 60, 90 seconds). Same as 31204 with LED and rechargeable power supply.



Supplied with head strip

Multilingual manual: GB, FR, IT, ES, PT, GR



**LED + BATTERY**

Suggested for short examinations

Adjustable illumination through a slight turn around of the switch

Multilingual manual: GB, FR, IT, ES, PT, GR



**LED + POWER SUPPLY**

Supplied with head strip, rechargeable power supply

Adjustable illumination

Manual: GB, IT

## USB VIDEONYSTAGMOSCOPE SYSTEM

- 31213 VIDEO-NYSTAGMOSCOPE - USB - mask with Ir camera and USB cable
- Facilitate to view and register the Nystagmus through standard vestibular tests (both position and instrumental proof). This system allows to view on the monitor or PC/Notebook images of patient eyes during tests and to make Visual Suppression Test (VST) through a green LED light inside the mask to allow fixation. Innovative technology allows to get quality images. Power 12V DC.

**MADE IN ITALY**



752x582 pixels resolution  
520 TV lines

31213