

## WAN STETHOSCOPES - ADULT



WAN STETHOSCOPES  
Chrome-plated

### WAN ADULT STETHOSCOPES

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring. Assure optimum acoustic performances.

Chestpiece: Ø 45 mm, bell Ø 30 mm.

Tubing: length 63 cm, int. Ø 4.2 mm, ext. Ø 8 mm.

Supplied in a cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.

GIMA code	WAN STETHOSCOPES - LATEX FREE	
32569	WAN single head adult - black	
32570	WAN single head adult - royal blue	
32571	WAN single head adult - burgundy	
32572	WAN single head adult - orange	
32574	WAN double head adult - black	
32575	WAN double head adult - royal blue	
32576	WAN double head adult - burgundy	
32577	WAN double head adult - orange	
32678	Diaphragm Ø 44 mm - spare	

## PEDIATRIC-NEONATAL STETHOSCOPES



### WAN PEDIATRIC AND NEONATAL STETHOSCOPES:

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring. Assure optimum acoustic performances. Tubing length: 63 cm.

Neonatal stethoscope: diaphragm Ø 22 mm, bell Ø 19 mm.

Pediatric stethoscope: diaphragm Ø 35 mm, bell Ø 30 mm.


Supplied in a cardboard box. Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.

### WAN PLUS:

Same as standard with inner spring binaural, superior tubing quality and 2 pair of soft colour coordinated ear tips in different size.

GIMA code	WAN STETHOSCOPES - LATEX FREE	
32508	WAN double head neonatal - red	
32509	WAN double head neonatal - grey	
32513	WAN double head pediatric - light-blue	
32515	WAN double head pediatric - red	
32511	WAN PLUS double head pediatric - pink	
32506	WAN PLUS double head pediatric - black	

## YTON STETHOSCOPES



### YTON STETHOSCOPES:


Gima Yton stethoscopes with single or double head.

Solid aluminium chestpiece Ø 45 mm, bell Ø 30 mm (for double head), with colour coordinated plastic retaining ring.

"Duracote" finished.

Available in four different colours: dark green, dark blue, red and aqua blue.

Manual and box in 11 languages: GB, FR, IT, DE, ES, PT, GR, PL, SE, HU, Arabic.



GIMA code	YTON STETHOSCOPES - LATEX FREE	
49500	Yton single head - dark green	
49501	Yton single head - dark blue	
49502	Yton single head - red	
49503	Yton single head - aqua blue	
49510	Yton double head - dark green	
49511	Yton double head - dark blue	
49512	Yton double head - red	
49513	Yton double head - aqua blue	